



EFFECTIVE
JUNE 26, 2017

1-800-MIDMARK
FAX 800-365-8631

Medical Product Guide

Midmark and Ritter Medical Products
Midmark Diagnostic Products

Product lines included:



Ritter[®]
by MIDMARK

General Information

Ordering Notes 4
Freight and Shipping..... 4
Midmark and Ritter Medical Equipment Warranty 5

Midmark Tables

Midmark 641 Barrier-Free® Power Procedures Chair 6-7
Midmark 630 HUMANFORM® Power Procedures Table 8-10
Midmark 647 Barrier-Free® Power Podiatry Procedures Chair..... 11-12
Midmark 646 Basic Power Podiatry Procedures Chair 13
Midmark 625 Barrier-Free® Power Examination Table 14-15
Midmark 604 Manual Examination Table..... 16-17
Midmark 640 Pediatric Examination Table..... 18
Midmark Upholstery Color Selections 19

Ritter Tables

Ritter 230 Power Procedures Table..... 20-21
Ritter 244 Barrier-Free® Bariatric Power Treatment Table..... 22
Ritter 223 Barrier-Free® Power Examination Table..... 23-24
Ritter 222 Barrier-Free® Power Examination Table..... 25-26
Ritter 204 Manual Examination Table 27-28
Ritter 203 Treatment Table 29
Ritter 95 Treatment Table 29
Ritter Upholstery Color Selections..... 30

Seating

Midmark 425 Air Lift Physician Stool (Hand Operated) 31
Midmark 427 Air Lift Physician Stool (Foot Operated) 31
Ritter 279 Air Lift Surgeon Stool (Hand Operated)..... 32
Ritter 278 Air Lift Surgeon Stool (Foot Operated) 32
Ritter 277 Air Lift Stool with Backrest..... 32
Ritter 276 Air Lift Stool..... 32
Ritter 275 Adjustable Physician Stool with Backrest..... 33
Ritter 274 Adjustable Physician Stool..... 33
Ritter 273 Air Lift Stool with Backrest..... 33
Ritter 272 Air Lift Stool..... 33
Ritter 271 Adjustable Stool with Backrest 34
Ritter 270 Adjustable Stool 34
Midmark 680 Upholstered Side Chair 35
Ritter 281 Blood Drawing Chair 35
Seating Upholstery Color Selections 36

Lighting

Ritter 255 LED Procedure Light 37
Ritter 253 LED Exam Light 38
Ritter 250 LED Exam Light 39

Diagnostic Products

Midmark Diagnostic Medical Products Warranty 40
Midmark Diagnostic Medical Products Support Information 41
Midmark Diagnostic Medical Products Ordering Information 42
Midmark IQvitals® and IQvitals® Zone™ and Supplies 43-47
Midmark IQspiro® Digital Spirometer, Software and Supplies 48
Midmark IQecg® Digital ECG, Software and Supplies..... 49
Midmark IQholter® Digital Holter, Software and Supplies 50-51
Midmark IQcart® Diagnostic Equipment Cart 52
Midmark Care Exchange® 6231 Procedure Workstation 52




Midmark Care Exchange® Workstations

Midmark Care Exchange® 6201 Traditional Style Laptop Workstation 53
 Midmark Care Exchange® 6202 Flat Panel PC Workstation 54
 Midmark Care Exchange® 6203 Flat Panel EMR Transition Workstation 55
 Midmark Care Exchange® 6204 Compact Rectangular Workstation 56
 Midmark Care Exchange® 6205 Compact Kidney Style Workstation..... 57
 Midmark Care Exchange® 6206 Kidney Style Laptop Workstation with Security System..... 58
 Midmark Care Exchange® 6207 Basic Kidney Style Laptop Workstation..... 59
 Midmark Care Exchange® 6208 Flat Panel Kidney Style Workstation 60
 Midmark Care Exchange® 6231 Procedure Workstation 61
 Midmark Care Exchange® 6261 AC Powered Traditional Style Laptop Workstation..... 62
 Midmark Care Exchange® 6263 AC Powered Flat Panel PC Workstation..... 63
 Midmark Care Exchange® 6281 Laptop Wall Mount 64
 Midmark Care Exchange® 6282 Standard Duty Wall Mount 64
 Midmark Care Exchange® 6283 Heavy Duty Wall Mount..... 65
 Midmark Care Exchange® 6271 Wall Mount Adj. Flip-Up Keyboard/Adj. Flat Screen..... 66
 Midmark Care Exchange® 6272 Wall Mount 67
 Midmark Care Exchange® 6273 Wall Mount 66
 Midmark Care Exchange® 6274 Wall Mount 68
 Midmark Care Exchange® 6275 Wall Mount 68
 Midmark Care Exchange® 6276 Wall Mount Adjustable Height 69



Synthesis® Casework and Treatment Cabinets by Midmark

Synthesis® Casework Warranty..... 70
 Synthesis® Casework Color Selections 71
 Synthesis® Exam Room Packages E1-E7 72-78
 Synthesis® Mobile Treatment Cabinets M2; M4-M5 79-81
 Synthesis® Glove/ Towel/ Cup Dispenser D1-D4..... 82
 Synthesis® Dressing Nook N1-N2..... 83
 Synthesis® Waste Receptacle W1..... 84
 Synthesis® Instrument Panel I1-I2 85
 Ritter 261 Waste Receptacle 86



Instrument Processing

Ritter M11 UltraClave® Automatic Sterilizer 87
 Ritter M9 UltraClave® Automatic Sterilizer 87
 Ritter M9D Manual Sterilizer..... 88
 Ritter M3 UltraFast® Automatic Sterilizer 88
 Sterilizer Accessories..... 89
 VistaCool™ Direct-to-Drain System..... 90
 QuickClean™ Ultrasonic Cleaner QC6/QC6R..... 91
 QuickClean™ Ultrasonic Cleaner QC3/QC3R..... 91
 QuickClean™ Ultrasonic Cleaner QC1 92
 Ultrasonic Parts & Accessories 92



Midmark Customer Service

To Place an Order – E-mail: Orders-Medical@midmark.com or Fax: 1-800-365-8631. For order status, requests and information, contact Medical Customer Service at:

1-844-856-1233

E-mail: allcustservice-medical@midmark.com.

Midmark Technical Service

For Service Part identification, technical/operational questions, warranty and product issues/complaints, contact Technical Service at: **1-844-856-1230.**

Sterilizer & Ultrasonic Products: sterilizers.ts@midmark.com

Diagnostic Products: techsupport@midmark.com

Procedure & Exam Tables, Seating & Lighting: examroomeq.ts@midmark.com

Medical Casework, Carts & Care Exchange Workstations: casework.ts@midmark.com

Special Notes

1. **ORDERING**—Please provide the following information on all orders:

- Purchase order number
- Bill to and ship to addresses
- Product model and/or part number(s)
- Description including upholstery color and/or panel/countertop color
- Quantity needed
- Pricing
- Required shipping date
- Special shipping requirements (see details below)
MINIMUM ORDER — \$100.

DROP SHIP REQUIREMENTS

- Please specify if the customer has a loading dock or has special shipping requirements.
- Receiving contact name, e-mail and phone number is required.

2. **SPECIAL DESIGN** — Orders for equipment of special design, construction, upholstery or special finish are not subject to cancellation after your acceptance of our acknowledgment.

3. **SPECIAL SHIPPING REQUIREMENTS** — These services are available to enhance delivery to your customers. Please call for pricing.

- **CALL BEFORE:** Carrier will call delivery contact to provide a drop-off timeframe, once they have possession of the product from Midmark and are at the final terminal.
- **LIFTGATE:** Carrier uses liftgate/ramp to get product to ground level when facility does not have a dock. Facility is responsible for taking product once at ground level.
- **INSIDE DELIVERY:** Carrier unloads product and delivers inside the first set of doors only, provided entrance to facility can accommodate freight as packaged (including skids) and there are no steps.
- **FULL SERVICE (WHITE GLOVE):** Carrier sets product in place per the delivery contact. Carrier is responsible for assembling product and removal of all packaging, trash, and skids from the delivery site (custom casework does not include installation).

4. **FREIGHT DAMAGE** — Merchandise arriving with external damage should be so described on bill of lading. In case of concealed damage, carrier must be notified within 15 days after receipt and claims for damage instituted.

5. **ORDER CHANGES** — Any changes made to an order within 48 hours of the scheduled ship date may affect the scheduled ship date, requiring Midmark to reschedule.

6. **ORDER CANCELLATION** — Any order, if cancelled within ten days before the scheduled shipping date, will be charged a 20% cancellation charge.

7. **RETURN OF MERCHANDISE** — Authorization must be obtained for the return of any item. If return is authorized, transportation charges must be prepaid. Risk of loss with respect to returned goods shall remain with purchaser until receipt by Midmark. Items returned for credit are subject to charge of 25% or more to cover inspection, damage in transit, reconditioning and repacking. Items may be returned within 6 months (except Service Parts within 90 days) from the original ship date if in a good, unused condition. RMA # should be written on the outside of the packaging. Diagnostic products may be returned within 30 days if in a good, unused condition and in original sealed packaging.

Note: Certain items are non-returnable: Specials, Synthesis® Casework, Care Exchange® Workstations, Discontinued Items or Parts, Custom or Non-standard Colors, Service Upholstery and Service Parts if value is less than \$50. Diagnostic Consumables e.g. Electrodes, Prep Kits, Pneumotach Mouthpieces, Temperature Probe Covers, Nose Clips, Thermal Paper, All Disposable Items and Diagnostic Items value less than \$50 per.

8. **SPECIAL ORDER REQUIREMENTS** — Any order requiring any special deviation from standard products will carry the following disclaimer:

We hereby accept and acknowledge the above captioned order, provided however, because of the special nature of this order in that this equipment is being supplied upon special design specifications, special materials and/or field modifications, as per your request, the standard warranties of Midmark and Ritter on its products are not applicable to this order, and this acceptance and acknowledgment is conditioned upon your indemnification and promise to hold Midmark and Ritter, its officers, employees and agents, harmless in respect to any and all warranty claims, damages, both direct and consequential, expenses, costs and any liability arising from the use of or relating to this special equipment.

9. **IMPROVE PRODUCT QUALITY** — Midmark Corporation reserves the right to change design and construction without notice.

Midmark and Ritter Medical Equipment (Tables, Physician Seating, Lights and Sterilizers) Limited Warranty

SCOPE OF WARRANTY Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will, at Midmark's option, repair or replace components of the domestic and international medical products manufactured by Midmark (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. The sole remedy under this limited warranty is the repair or replacement, at Midmark's option, of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which are determined to exist upon examination by Midmark. This warranty extends only to the original retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

Midmark warrants to the original retail purchaser that during the applicable warranty period it will repair or replace software contained within the products manufactured by Midmark (except for those not warranted under "Exclusions") if: (1) the media on which the software is furnished exhibits defects in material or workmanship under normal use; or (2) the software does not substantially conform to its published specifications.

APPLICABLE WARRANTY PERIOD The applicable warranty period, measured from the date of invoice to the original retail purchaser of the product and shall be one (1) year for all warranted products and components.

OBTAINING WARRANTY SERVICE Warranty service must be obtained through either Midmark or an authorized dealer in the Midmark product line for which warranty service is requested. Midmark may be contacted for warranty service inquiries or issues via email at www.midmark.com; by phone at 1-800-MIDMARK or by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380.

It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this limited warranty.

EXCLUSIONS: This limited warranty does not cover and Midmark shall not be liable for the following:

- (1) Defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) Products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environmental conditions and electrical requirements;
- (3) Products considered to be of a consumable or sterile nature;
- (4) Accessories or parts not manufactured by Midmark;
- (5) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (6) Costs and expenses of routine maintenance and cleaning; and
- (7) Representations and warranties made by any person or entity other than Midmark.
- (8) Matching of color, grain or texture except to commercially acceptable standards;
- (9) Changes in color caused by natural or artificial light;
- (10) Custom manufactured products;
- (11) Alterations or modifications to the product by any person or entity other than Midmark; and
- (12) Products that would otherwise be covered under this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical market.

SOFTWARE; WITH RESPECT TO SOFTWARE THAT IS A PRODUCT OR COMPONENT THEREOF, MIDMARK DOES NOT WARRANT THAT THE SOFTWARE: (1) IS ERROR FREE; (2) CAN BE USED WITHOUT PROBLEMS OR INTERRUPTIONS; OR (3) IS FREE FROM VULNERABILITY TO INTRUSION OR ATTACK BY VIRUSES OR OTHER METHODS.

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, LOSS OF DATA, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. THIS DISCLAIMER SHALL SURVIVE ANY FAILURE OR ASSERTED FAILURE OF THE ESSENTIAL PURPOSE OF THIS LIMITED WARRANTY OR ITS REMEDIES SPECIFIED HEREIN.

WARRANTY DISCLAIMER THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

STATUTE OF LIMITATIONS No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

NO AUTHORIZATION No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

*Warranty information included in product manuals



Midmark 641

Barrier-Free® Power Procedures Chair

- Power height, back, foot and tilt
- 450-pound patient weight capacity
- 19" to 40" seat height adjustment (non-rotation models)
- 22.5" to 43.5" seat height adjustment (rotation models)
- Perimeter contouring for maximum access to patient
- Flip-out foot extension
- Hand and foot control standard
- Paper roll holder standard
- Choice of Premium seamless or Ultra-Premium plush upholstery
- Meets UL 60601-1 and CAN/CSA C22.2 No. 601.1-M90 standards
- Table weight: 450 lb
- Shipping weight: 500 lb

Model	Description
Standard Models	
641-002 002-0966-00-XXX ¹	Procedures Chair Non-Programmable w/o receptacle, w/o top Upholstery top - Premium
641-002 002-0968-00-XXX ¹	Procedures Chair Non-Programmable w/o receptacle, w/o top Upholstery top - Ultra-Premium
641-003 002-0966-00-XXX ¹	Procedures Chair Programmable w/receptacle, w/o top Upholstery top - Premium
641-003 002-0968-00-XXX ¹	Procedures Chair Programmable w/receptacle, w/o top Upholstery top - Ultra-Premium
Rotation Models	
641-004 002-0966-00-XXX ¹	Procedures Chair Non-Programmable w/rotation, w/o top Upholstery top - Premium
641-004 002-0968-00-XXX ¹	Procedures Chair Non-Programmable w/rotation, w/o top Upholstery top - Ultra-Premium
641-005 002-0966-00-XXX ¹	Procedures Chair Programmable w/rotation, w/o top Upholstery top - Premium
641-005 002-0968-00-XXX ¹	Procedure Chair Programmable w/rotation, w/o top Upholstery top - Ultra-Premium

Accessories	Description
Arms/Armboards	
9A179002-XXX ^{1,6}	Fixed armboard
9A396001-XXX ^{1,4,6}	Chair arms, Premium upholstery
9A398001-XXX ^{1,4,6}	Chair arms, Ultra-Premium upholstery
9A81001-XXX ^{1,6}	Articulating armboard
9A82001-XXX ^{1,6}	Hand surgery armboard
Casters	
9A330001 ²	Caster base, swivel wheels
9A365001 ²	Casters
Controls	
9A408001 ³	Wireless control, non-programmable, linear-style, plant installed
9A408002 ³	Wireless control, programmable, linear-style, plant installed
9A408003 ³	Wireless control, non-programmable, linear-style, field installed
9A408004 ³	Wireless control, programmable, linear-style, field installed
Headrests	
9A378001-XXX ^{1,6}	Flat headrest, Premium upholstery
9A380001-XXX ^{1,6}	Flat headrest, Ultra-Premium upholstery
9A382001-XXX ^{1,6}	U-shaped headrest, Premium upholstery
9A384001-XXX ^{1,6}	U-shaped headrest, Ultra-Premium upholstery
9A386001-XXX ^{1,6}	Oval-shaped headrest, Premium upholstery
9A388001-XXX ^{1,6}	Oval-shaped headrest, Ultra-Premium upholstery
9A390001-XXX ^{1,6}	Rectangular-shaped headrest, Premium upholstery
9A392001-XXX ^{1,6}	Rectangular-shaped headrest, Ultra-Premium upholstery
9A394001-XXX ^{1,6}	Magnetic headrest, Premium upholstery
9A395001-XXX ^{1,6}	Magnetic headrest, Ultra-Premium upholstery
Instrument Trays/Work Shelves/IV Poles	
9A427001	Swing arm instrument tray, fixed top
9A427002	Swing arm instrument tray, extension top
9A428001	Swing arm PC workshelf
9A267003 ⁵	Instrument tray, double arm
9A77001 ⁵	IV pole
Miscellaneous	
9A152010	3' power cord
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A406001	Security straps, plant installed
9A406002	Security straps, field installed
9A78001-XXX ^{1,5,6}	Vision block screen
9A416001	Foot section plastic cover
002-0966-00-999	Upholstery top - Premium Special Color
002-0968-00-999	Upholstery top - Ultra-Premium Special Color
002-0967-00-XXX ¹	CAL 133 Upholstery top - Premium
002-0969-00-XXX ¹	CAL 133 Upholstery top - Ultra-Premium
002-0967-00-999	CAL 133 Upholstery top - Premium Special Color
002-0969-00-999	CAL 133 Upholstery top - Ultra-Premium Special Color

¹ - See color code chart on page 19

² - Not available on rotation models

³ - Accessory supplied with hand and foot control (batteries included)

⁴ - Interference with armboard accessories will likely occur

⁵ - Requires 9A20500X base rail

⁶ - Additional charge for Special or Ultra Premium Upholstery Color for Accessory



Midmark 630

HUMANFORM® Power Procedures Table

- Power height, back, foot and tilt
- Motion Profiling® designed to provide more natural, human-like movement that works with the patient and provider for a clear view of the surgical site and improved comfort
- Active Sensing Technology® (AST) is designed to actively protect patient and provider by sensing and avoiding harmful impact during movement, such as objects being caught under the seat and footrest
- Barrier-Free® design
- 650-pound patient weight capacity
- 18" to 40" seat height adjustment (non-rotational models)
- 21.5" to 43.5" seat height adjustment (rotational models)
- Hand and foot controls standard
- Available in 28" or 32" upholstery widths
- Meets UL 60601-1 and CAN/CSA C22.2 No.601.1-M90 standards
- Table weight: 450 lb
- Shipping weight: 525 lb

Model	Description
Non-Programmable	
630-010 002-1644-XXX ¹	Procedures Table Non-Programmable w/receptacle, w/o top Upholstery top - 28" wide Premium
630-010 002-1646-XXX ¹	Procedures Table Non-Programmable w/receptacle, w/o top Upholstery top - 28" wide Ultra-Premium
630-011 002-1644-XXX ¹	Procedures Table Non-Programmable wireless, w/receptacle, w/o top Upholstery top - 28" wide Premium
630-011 002-1646-XXX ¹	Procedures Table Non-Programmable wireless, w/receptacle, w/o top Upholstery top - 28" wide Ultra-Premium
Non-Programmable with Rotation	
630-012 002-1644-XXX ¹	Procedures Table Non-Programmable w/rotation, w/o top Upholstery top - 28" wide Premium
630-012 002-1646-XXX ¹	Procedures Table Non-Programmable w/rotation, w/o top Upholstery top - 28" wide Ultra-Premium
630-013 002-1644-XXX ¹	Procedures Table Non-Programmable wireless, w/rotation, w/o top Upholstery top - 28" wide Premium
630-013 002-1646-XXX ¹	Procedures Table Non-Programmable wireless, w/rotation, w/o top Upholstery top - 28" wide Ultra-Premium

Model	Description
Programmable 630-020 002-1644-XXX ¹	Procedures Table Programmable w/receptacle, w/o top Upholstery top - 28" wide Premium
630-020 002-1646-XXX ¹	Procedures Table Programmable, w/receptacle, w/o top Upholstery top - 28" wide Ultra-Premium
630-021 002-1644-XXX ¹	Procedures Table Programmable wireless, w/receptacle, w/o top Upholstery top - 28" wide Premium
630-021 002-1646-XXX ¹	Procedures Table Programmable wireless, w/receptacle, w/o top Upholstery top - 28" wide Ultra-Premium
Programmable with Rotation	
630-022 002-1644-XXX ¹	Procedures Table Programmable, w/rotation, w/o top Upholstery top - 28" wide Premium
630-022 002-1646-XXX ¹	Procedures Table Programmable, w/rotation, w/o top Upholstery top - 28" wide Ultra-Premium
630-023 002-1644-XXX ¹	Procedures Table Programmable wireless, w/rotation, w/o top Upholstery top - 28" wide Premium
630-023 002-1646-XXX ¹	Procedures Table Programmable wireless, w/rotation, w/o top Upholstery top - 28" wide Ultra-Premium
Upholstery Options	
002-1644-999 002-1646-999	Upholstery top - 28" wide Premium Special Color Upholstery top - 28" wide Ultra-Premium Special Color
002-1648-XXX ¹ 002-1648-999	Upholstery top - 32" wide Premium Upholstery top - 32" wide Premium Special Color
002-1650-XXX ¹ 002-1650-999	Upholstery top - 32" wide Ultra-Premium Upholstery top - 32" wide Ultra-Premium Special Color
002-1645-XXX ¹ 002-1645-999	CAL 133 Upholstery top - 28" wide Premium CAL 133 Upholstery top - 28" wide Premium Special Color
002-1647-XXX ¹ 002-1647-999	CAL 133 Upholstery top - 28" wide Ultra-Premium CAL 133 Upholstery top - 28" wide Ultra-Premium Special Color
002-1649-XXX ¹ 002-1649-999	CAL 133 Upholstery top - 32" wide Premium CAL 133 Upholstery top - 32" wide Premium Special Color
002-1651-XXX ¹ 002-1651-999	CAL 133 Upholstery top - 32" wide Ultra-Premium CAL 133 Upholstery top - 32" wide Ultra-Premium Special Color

Midmark Product Guide

Effective June 26, 2017

Accessories	Description
Arm Systems/Armboards	
9A550001-XXX ^{1,6}	Chair arms, 28" wide Premium upholstery
9A550003-XXX ^{1,6}	Chair arms, 28" wide Ultra-Premium upholstery
9A550005-XXX ^{1,6}	Chair arms, 32" wide Premium upholstery
9A550007-XXX ^{1,6}	Chair arms, 32" wide Ultra-Premium upholstery
9A179002-XXX ^{1,6,3}	Fixed Armboard
9A81001-XXX ^{1,6,3}	Articulating Armboard
9A82001-XXX ^{1,6,3}	Hand surgery Armboard
Controls (additional)	
9A552001	Wired hand control, non-programmable, field installed
9A553001	Wired foot control, non-programmable, field installed
9A552002	Wired hand control, programmable, field installed
9A553002	Wired foot control, programmable, field installed
Headrests	
9A496001-XXX ^{1,6}	Rectangular-shaped headrest, Premium upholstery
9A497001-XXX ^{1,6}	Rectangular-shaped headrest, Ultra-Premium upholstery
9A500001-XXX ^{1,6}	Oval-shaped headrest, Premium upholstery
9A501001-XXX ^{1,6}	Oval-shaped headrest, Ultra-Premium upholstery
9A504001-XXX ^{1,6}	Round headrest, Premium upholstery
9A505001-XXX ^{1,6}	Round headrest, Ultra-Premium upholstery
9A508001-XXX ^{1,6}	Flat headrest, Premium upholstery
9A509001-XXX ^{1,6}	Flat headrest, Ultra-Premium upholstery
9A512001-XXX ^{1,6}	U-shaped headrest, Premium upholstery
9A513001-XXX ^{1,6}	U-shaped headrest, Ultra-Premium upholstery
9A516001-XXX ^{1,6}	Magnetic headrest, Premium upholstery
9A517001-XXX ^{1,6}	Magnetic headrest, Ultra-Premium upholstery
Instrument Trays/Work Shelves/IV Poles	
9A363001	PC Work Surface
9A267003 ³	Instrument tray, double arm
9A77001 ³	IV pole
Knee Crutches	
9A410001 ⁴	Standard knee crutches (crutches only)
9A410002 ⁴	Standard knee crutch system
9A411001 ⁴	Articulating knee crutches (crutches only)
9A411002 ⁴	Articulating knee crutch system
9A412001 ⁴	Knee crutch mounting brackets (brackets only)
Accessory Rails	
9A521001	Seat Accessory Rails, Surgical Size, 28" wide upholstery, plant installed
9A522001	Seat Accessory Rails, Surgical Size, 32" wide upholstery, plant installed
9A567001	Seat Accessory Rails, Surgical Size, 28" wide upholstery, field installed
9A568001	Seat Accessory Rails, Surgical Size, 32" wide upholstery, field installed
9A557001	Accessory Rails, 28" wide upholstery, field installed
9A558001	Accessory Rails, Surgical Size, 28" wide upholstery, field installed
9A557002	Accessory Rails, 32" wide upholstery, field installed
9A558002	Accessory Rails, Surgical Size, 32" wide upholstery, field installed
Miscellaneous	
9A554001 ²	Casters - set of 4
9A152010	3' power cord
9A157003 ³	Security side panels
9A559001 ⁵	Security straps, field installed
9A555001	Gynecology debris tray accessory
9A556001	Gynecology debris tray (tray only)
9A78001-XXX ^{1,3,6}	Anesthesia / Vision block screen
AD 1779 ²	Tamper resistant duplex receptacle

¹ - See color code chart on page 19

² - Not available on rotation models

³ - Requires (1) 9A557001/9A557002/9A558001/9A558002 accessory rails

⁴ - Not compatible with 9A521001, 9A522001, 9A567001, or 9A568001

⁵ - Requires (1) 9A557001/9A557002/9A558001/9A558002 accessory rails and (1) 9A521001/9A522001 seat rails

⁶ - Additional charge for Special or Ultra Premium Upholstery Color for Accessory

Midmark 647

Barrier-Free® Power Podiatry Procedures Chair



- Power height, tilt and back
- 450-pound patient weight capacity
- 19" to 40" seat height adjustment (non-rotation models)
- 22.5" to 43.5" seat height adjustment (rotation models)
- 9" foot section extension
- Sealed foot control standard
- Integrated touch controls standard
- Chair arms standard
- Choice of Premium seamless or Ultra-Premium plush upholstery
- Protective foot cover standard with upholstery
- Meets UL 60601-1 and CAN/CSA C22.2 No. 601.1-M90 standards
- Table weight: 433 lb (non-rotation models)
- Shipping weight: 515 lb (non-rotation models)
- Table weight: 539 lb (rotation models)
- Shipping weight: 621 lb (rotation models)

Model	Description
647-001 002-1119-XXX ¹	Barrier-Free® Podiatry Chair Non-programmable w/receptacles, w/o top Upholstery top - Premium
647-001 002-1121-XXX ¹	Barrier-Free® Podiatry Chair Non-programmable w/receptacles, w/o top Upholstery top - Ultra-Premium
647-002 002-1119-XXX ¹	Barrier-Free® Podiatry Chair Programmable w/receptacles, w/o top Upholstery top - Premium
647-002 002-1121-XXX ¹	Barrier-Free® Podiatry Chair Programmable w/receptacles, w/o top Upholstery top - Ultra-Premium
647-003 002-1119-XXX ¹	Barrier-Free® Podiatry Chair Non-programmable w/rotation, w/o receptacles & top Upholstery top - Premium
647-003 002-1121-XXX ¹	Barrier-Free® Podiatry Chair Non-programmable w/rotation, w/o receptacles & top Upholstery top - Ultra-Premium
647-004 002-1119-XXX ¹	Barrier-Free® Podiatry Chair Programmable w/rotation, w/o receptacles & top Upholstery top - Premium
647-004 002-1121-XXX ¹	Barrier-Free® Podiatry Chair Programmable w/rotation, w/o receptacles & top Upholstery top - Ultra-Premium

Accessories	Description
Instrument Trays/Work Shelves/IV Poles	
9A427001	Swing arm instrument tray, fixed top
9A427002	Swing arm instrument tray, extension top
9A428001	Swing arm PC work shelf
9A432001	Instrument tray, swing arm (tray only)
9A77001 ³	IV Pole
Casters	
9A330001 ²	Caster base, swivel wheels
9A365001 ²	Casters
Miscellaneous	
9A418001	Wireless foot control, non-programmable, linear-style, plant installed
9A418002	Wireless foot control, non-programmable, linear-style, field installed
9A418003	Wireless foot control, programmable, linear-style, plant installed
9A418004	Wireless foot control, programmable, linear-style, field installed
9A422001-XXX ^{1,4}	Vision block screen
9A420001	Plastic protective foot cover
9A423001	Paper roll holder
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A152010	3' power cord
002-1119-999	Upholstery top - Premium Special Color
002-1121-999	Upholstery top - Ultra-Premium Special Color
002-1120-XXX ¹	CAL 133 Upholstery top - Premium
002-1122-XXX ¹	CAL 133 Upholstery top - Ultra-Premium
002-1120-999	CAL 133 Upholstery top - Premium Special Color
002-1122-999	CAL 133 Upholstery top - Ultra-Premium Special Color

¹ - See color code chart on page 19

² - Not available on rotation models

³ - Requires 9A20500X base rail

⁴ - Additional charge for Special or Ultra Premium Upholstery Color for Accessory



Midmark 646

Basic Power Podiatry Procedures Chair

- Power tilt and back
- 450-pound patient weight capacity
- 21" fixed seat height
- 9" foot section extension
- Sealed foot control standard
- Integrated touch controls standard
- Chair arms standard
- Choice of Premium seamless or Ultra-Premium plush upholstery
- Protective foot cover standard with upholstery
- Meets UL 60601-1 and CAN/CSA C22.2 No. 601.1-M90 standards
- Table weight: 323 lb
- Shipping weight: 405 lb

Model	Description
646-001 002-1119-XXX ¹	Podiatry Chair Non-programmable w/receptacles, w/o top Upholstery top - Premium
646-001 002-1121-XXX ¹	Podiatry Chair Non-programmable w/receptacles, w/o top Upholstery top - Ultra-Premium

Accessories	Description
Instrument Trays/Work Shelves/IV Poles	
9A427001	Swing arm instrument tray, fixed top
9A427002	Swing arm instrument tray, extension top
9A428001	Swing arm PC work shelf
9A432001	Instrument tray, swing arm (tray only)
9A77001 ²	IV Pole
Casters	
9A330001	Caster base, swivel wheels
9A365001	Casters
Miscellaneous	
9A422001-XXX ^{1,3}	Vision block screen
9A420001	Plastic protective foot cover
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A152010	3' power cord
002-1119-999	Upholstery top - Premium Special Color
002-1121-999	Upholstery top - Ultra-Premium Special Color
002-1120-XXX ¹	CAL 133 Upholstery top - Premium
002-1122-XXX ¹	CAL 133 Upholstery top - Ultra-Premium
002-1120-999	CAL 133 Upholstery top - Premium Special Color
002-1122-999	CAL 133 Upholstery top - Ultra-Premium Special Color

¹ - See color code chart on page 19

² - Requires 9A20500X base rail

³ - Additional charge for Special or Ultra Premium Upholstery Color for Accessory



Midmark 625

Barrier-Free® Power Examination Table

- Power height and back
- 650-pound patient weight capacity
- 18" to 37" seat height adjustment
- Unique Quick Exam® and Motion Profiling® for efficient positioning standard
- Hand and foot control standard
- LH and RH receptacles
- Drawer heater (35 watt) - approx. temperature of 98°F (37°C)
- Choice of Premium seamless or Ultra-Premium plush upholstery
- Meets UL 60601-1 and CAN/CSA C22.2 No. 601.1-M90 standards
- Table weight: 433 lb
- Shipping weight: 505 lb

Model	Description
Standard Models	
625-001	Barrier-Free® Examination Table w/receptacles, drawer heater, pelvic tilt, w/o top
002-0874-XXX ¹	Upholstery top - 28" wide Premium
625-001	Barrier-Free® Examination Table w/receptacles, drawer heater, pelvic tilt, w/o top
002-1180-XXX ¹	Upholstery top - 32" wide Premium
625-001	Barrier-Free® Examination Table w/receptacles, drawer heater, pelvic tilt, w/o top
002-0875-XXX ¹	Upholstery top - 28" wide Ultra-Premium
625-003	Barrier-Free® Examination Table w/receptacles, drawer heater, pelvic tilt, w/wireless controls, w/o top
002-0874-XXX ¹	Upholstery top - 28" wide Premium
625-003	Barrier-Free® Examination Table w/receptacles, drawer heater, pelvic tilt, w/wireless controls, w/o top
002-1180-XXX ¹	Upholstery top - 32" wide Premium
625-003	Barrier-Free® Examination Table w/receptacles, drawer heater, pelvic tilt, w/wireless controls, w/o top
002-0875-XXX ¹	Upholstery top - 28" wide Ultra-Premium
625-004	Barrier-Free® Examination Table w/receptacles, drawer heater, pelvic tilt, digital scale, USB hub, wireless controls, w/o top
002-0874-XXX ¹	Upholstery top - 28" wide Premium
625-004	Barrier-Free® Examination Table w/receptacles, drawer heater, pelvic tilt, digital scale, USB hub, wireless controls, w/o top
002-1180-XXX ¹	Upholstery top - 32" wide Premium
625-004	Barrier-Free® Examination Table w/receptacles, drawer heater, pelvic tilt, digital scale, USB hub, wireless controls, w/o top
002-0875-XXX ¹	Upholstery top - 28" wide Ultra-Premium

Accessories	Description
Arm Systems	
9A561001 ^{4,6}	Patient Support Rails, rails only, pair, field installed
9A585001 ⁷	Patient Support Rail Receivers - pair, plant installed
9A582001 ⁵	Patient Support Rail Receiver Extensions - pair, field installed, for 32" wide upholstery
9A314001 ²	Adjustable arm system - arms only, field installed
9A314004 ⁸	Adjustable arm system - accessory receiver only 28" wide upholstery, plant installed
9A314005	Adjustable arm system - accessory receiver only 28" wide upholstery, field installed
9A314006 ⁸	Adjustable arm system - accessory receiver only 32" wide upholstery, plant installed
9A314007	Adjustable arm system - accessory receiver only 32" wide upholstery, field installed
9A354001 ³	Assist handles (handles only) 28" wide upholstery only, field installed
9A354002	Assist handles receiver only, 28" wide upholstery only, plant installed
9A354003	Assist handles receiver only, 28" wide upholstery only, field installed
Knee Crutches	
9A410001	Standard knee crutch only (204/604/22X/62X/230/630) field installed
9A410006	Standard knee crutch kit with accessory receiver (625) - field installed
9A411001	Articulating knee crutch only (204/604/22X/62X/230/630) field installed
9A411006	Articulating knee crutch kit with accessory receiver (625) - field installed
9A412005	Mounting brackets, hardware with accessory receiver (625) - field installed
Vital Signs Accessories	
9A478001	Midmark IQscale [®] serial cable for IQvitals [®] , 6 foot length, coiled
9A478002	Midmark IQscale [®] serial cable for IQvitals [®] , 15 foot length, straight
9A478003	Midmark IQscale [®] serial cable for IQvitals [®] , 30 foot length, straight
9A478004	Midmark IQscale [®] serial cable for IQvitals [®] , 50 foot length, straight
9A520001	USB A-B cable, 15 foot length
9A551001	IQvitals [®] Table Mount, IQhub [®]
9A551002	IQvitals [®] Table Mount, Non-IQhub [®]
9A479001	Welch Allyn [®] Spot Vitals LXi serial cable, 6 foot length, coiled
9A479002	Welch Allyn [®] Spot Vitals LXi serial cable, 15 foot length, straight
002-10063-00	Connectivity Wiring Kit (for 625-004)
Miscellaneous	
9A104001	Urology drain pan with hose
9A363001	PC Work Surface
9A429010	Wireless Control
9A431001	Supply assistant storage unit
9A70000	Stainless steel treatment pan
002-0874-999	Upholstery top - 28" Wide Premium Special Color
002-1180-999	Upholstery top - 32" Wide Premium Special Color
002-0875-999	Upholstery top - 28" Wide Ultra-Premium Special Color
002-0878-XXX ¹	CAL 133 Upholstery Top - 28" Wide Premium
002-0879-XXX ¹	CAL 133 Upholstery Top - 28" Wide Ultra-Premium
002-1181-XXX ¹	CAL 133 Upholstery Top - 32" Wide Premium
002-0878-999	CAL 133 Upholstery Top - 28" Wide Premium Special Color
002-0879-999	CAL 133 Upholstery Top - 28" Wide Ultra-Premium Special Color
002-1181-999	CAL 133 Upholstery Top - 32" Wide Premium Special Color

¹ - See color code chart on page 19

² - 9A314004, 9A314005 or 9A314007 arm receiver required

³ - 9A354003 or 9A354002 assist handle receiver required

⁴ - 9A585001 patient support rail receivers required

⁵ - 9A582001 receiver extensions required for 32" wide upholstery

⁶ - If also mounting exam light to table, use 253-003 to minimize contact point

⁷ - Cannot order with 9A314004, 9A314006 or 9A354002

⁸ - Cannot order with 9A585001



Midmark 604

Manual Examination Table

- 500-pound patient weight capacity
- Electrical receptacle standard
- Choice of Premium seamless or Ultra-Premium plush upholstery
- Unique flex storage design, ships as RH, easily field configured to LH
- Manual back with air spring assist
- Meets UL 60601-1 and CAN/CSA C22.2 No.601.1-M90 standards
- Please specify upholstery/insert color
- Table weight: 250 lb
- Shipping weight: 300 lb

Model	Description
Standard Models	
604-001 002-0874-XXX ¹	Manual Examination Table w/o top Upholstery top - Premium
604-001 002-0875-XXX ¹	Manual Examination Table w/o top Upholstery top - Ultra-Premium
604-002 002-0874-XXX ¹	Manual Examination Table w/drawer heater, pelvic tilt w/o top Upholstery top - Premium
604-002 002-0875-XXX ¹	Manual Examination Table w/drawer heater, pelvic tilt w/o top Upholstery top - Ultra-Premium
Heated Top Models	
604-003 002-0876-XXX ¹	Manual Examination Table w/drawer heater, heated upholstery prepared pelvic tilt w/o top Heated upholstery top - Premium
604-003 002-0877-XXX ¹	Manual Examination Table w/drawer heater, heated upholstery prepared pelvic tilt w/o top Heated upholstery top - Ultra-Premium

Accessories	Description
Arm/Handle Systems	
9A314001 ²	Adjustable arm system – arms only, field installed
9A354001 ²	Assist handles (handles only), 28" wide upholstery only, field installed
9A354002	Accessory receiver only, 28" wide upholstery only, plant installed
9A354003	Accessory receiver only, 28" wide upholstery only, field installed
Knee Crutches³	
9A410001	Standard Knee Crutch only (204/604/22X/62X/230/630) - field installed
9A410004	Standard Knee Crutch Kit with accessory receiver (204/604/22X/62X) - field installed
9A411001	Articulating Knee Crutch only (204/604/22X/62X/230/630) - field installed
9A411004	Articulating Knee Crutch Kit with accessory receiver (204/604/22X/62X) - field installed
9A412003	Mounting Brackets, hardware with accessory receiver (204/604/22X/62X) - field installed
9A413003-271	604 cladding only for knee crutches, Fossil, field installed
Miscellaneous	
9A104001	Urology drain pan w/hose
9A363001	PC Work Surface Kit
9A70000	Stainless steel treatment pan
9A551002	IQvitals® Table Mount, Non-IQhub®
002-0874-999	Upholstery Top - Premium Special Color
002-0875-999	Upholstery Top - Ultra-Premium Special Color
002-0876-999	Heated Upholstery Top - Premium Special Color
002-0877-999	Heated Upholstery Top - Ultra-Premium Special Color
002-0878-XXX ¹	CAL 133 Upholstery Top - 28" Wide Premium
002-0879-XXX ¹	CAL 133 Upholstery Top - 28" Wide Ultra-Premium
002-0878-999	CAL 133 Upholstery Top - 28" Wide Premium Special Color
002-0879-999	CAL 133 Upholstery Top - 28" Wide Ultra-Premium Special Color

¹ - See color code chart on page 19

² - 9A354002 or 9A354003 accessory receiver required

³ - 9A413003-271 cladding required



Midmark 640 Pediatric Examination Table


- Pebble grey base and panel color only
- An upholstered end panel and beveled sides designed to help keep infants from rolling off
- Flat sides with no sharp edges
- Accurate and easy to use, built-in scale
- 36-inch exam surface height
- Ample storage and a hidden paper roll holder
- Meets ANSI/AAMI ES 60601-1:2005 & CAN/CSA C22.2 No. 60601-1-08 standards
- 42" W x 51.5" H x 23"D
- Table weight: 170 lb
- Shipping weight: 210 lb

Model	Description
640-001-XXX ^{1,2}	Pediatric Examination Table w/digital scale
640-001-999 ²	Pediatric Examination Table w/digital scale - Special Color
Miscellaneous	
640-011-XXX ^{1,2}	CAL 133 Pediatric Examination Table w/digital scale
640-012-XXX ^{1,2}	CAL 133 Pediatric Examination Table with door and drawer lock, plant installed
640-011-999 ²	CAL 133 Pediatric Examination Table w/digital scale - Special Color
640-012-999 ²	CAL 133 Pediatric Examination Table w/door and drawer lock, plant installed - Special Color
640-002-XXX ^{1,2}	Pediatric Examination Table w/door and drawer lock, plant installed
640-002-999 ²	Pediatric Examination Table w/door and drawer lock, plant installed - Special Color

¹ - See color code chart on page 19

² - Ultra-Premium colors are not available on the Midmark 640

Midmark Color-Specific Product Number Chart



Example:

630-010 *Table*

002-1644-235 *Top*

630 HUMANFORM® Power Procedures Table - Blueberry, Premium Upholstery Top

604-001 *Table*

002-0875-245 *Top*

604 Full-Featured Examination Table - Ultra-Sky, Ultra-Premium Upholstery Top & Inserts

Solid Premium Colors

235	Blueberry
236	Silver Sage
237	Tea Green
238	Fossil Grey
239	Cashmere Blue
240	Terra Cotta
312	Black

Pattern Premium Colors

247	Belagio
248	Siena
249	Tuscany
250	Firenze
251	Capri
252	Treviso

Ultra-Premium Colors¹

228	Ultra-Raven Wing
241	Ultra-Pacific Blue
242	Ultra-Spruce
243	Ultra-Olivine
244	Ultra-Arctic
245	Ultra-Sky
246	Ultra-Canyon

¹ See product price sheet for Ultra-Premium color availability and price.



Ritter 230

Power Procedures Table

- Power height, back, foot and tilt
- 450-pound patient weight capacity
- 22.5" to 40" seat height adjustment (non-rotation models)
- 26" to 43.5" seat height adjustment (rotation models)
- Knee-chest configurable
- Hand control standard
- Meets UL60601-1 and CAN/CSA C22.2 No. 601.1 – M90 standards
- Table weight: 455 lb
- Shipping weight: 525 lb

Model	Description
230-001 002-0824-XXX ¹	Procedures Table w/o top Seamless Upholstery top
230-002 002-0824-XXX ¹	Procedures Table w/receptacle, w/o top Seamless Upholstery top
230-003 002-0824-XXX ¹	Procedures Table w/rotation, w/o top Seamless Upholstery top

¹ - See color code chart on page 30

Accessories	Description
Arm Systems/Armboards	
9A179002-XXX ^{1,5}	Fixed armboard
9A328001	Adjustable arm system, plant installed
9A328002	Adjustable arm system, field installed
9A350001	Adjustable arm system (receiver only), plant installed
9A350002	Adjustable arm system (receiver only), field installed
9A350003	Adjustable arm system (arms only)
9A81001-XXX ^{1,5}	Articulating armboard
9A82001-XXX ^{1,5}	Hand surgery armboard
Casters	
9A330001 ^{2,3}	Caster base, swivel wheels
9A365001 ³	Casters
Controls	
9A371003	Foot control, linear style
Headrests/Footrests	
9A337001-XXX ^{1,5}	Larger 24" wide footrest
9A51001-XXX ^{1,5}	Facial pad
Instrument Trays/Work Shelves/IV Poles	
9A427001	Swing arm instrument tray, fixed top
9A427002	Swing arm instrument tray, extension top
9A428001	Swing arm PC workshelf
9A267003 ⁴	Instrument tray, double arm
9A77001 ⁴	IV pole
Knee Crutches	
9A410001	Standard knee crutches (crutches only)
9A410002	Standard knee crutch system
9A411001	Articulating knee crutches (crutches only)
9A411002	Articulating knee crutch system
9A412001	Knee crutch mounting brackets (brackets only)
Miscellaneous	
002-0950-XXX ¹	Flat upholstery set, 24" width
9A152010	3' Power Cord
9A157003	Security side panels
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A214001	Security straps, plant installed
9A214002	Security straps, field installed
9A329003 ³	Urology drain system, plant installed
9A329004 ³	Urology drain system, field installed
9A78001-XXX ^{1,4,5}	Vision block screen
002-0824-999	Seamless Upholstery top Special Color
002-0825-XXX ¹	CAL 133 Upholstery Top-Premium
002-0951-XXX ¹	CAL 133 Flat Upholstery set, 24" width
002-0825-999	CAL 133 Upholstery Top - Seamless Special Color
002-0951-999	CAL 133 Flat Upholstery set, 24" width Special Color
002-0950-999	Flat Upholstery set, 24" width Special Color

¹ - See color code chart on page 30

² - Not compatible with 9A337001

³ - Not available on rotation models

⁴ - Requires 9A20500X base rail

⁵ - Additional charge for Special or Premium Upholstery Color for Accessory



Ritter 244

Barrier-Free® Bariatric Power Treatment Table

- Power height, back and tilt
- Reverse trendelenburg
- 850-pound patient weight capacity
- 18" to 34" seat height adjustment
- 32" wide seat
- Hand control standard
- Meets UL 60601-1 and CAN/CSA C22.2 No.601.1-M90 standards
- Table weight: 340 lb
- Shipping weight: 400 lb

Model	Description
244-001 002-0861-XXX ¹	Barrier-Free® Bariatric Table w/o top Seamless Upholstery top
Accessories	
9A314001	Adjustable arm system, field installed
9A346001	Casters
9A347001	Foot control
9A152010	3' Power Cord
Miscellaneous	
002-0861-999	Seamless Upholstery Top Special Color
002-1017-XXX ¹	CAL 133 Seamless Upholstery Top
002-1017-999	CAL 133 Seamless Upholstery Top Special Color

¹ - See color code chart on page 30



Ritter 223

Barrier-Free® Power Examination Table

- Power height and back
- 400-pound patient weight capacity
- 18" to 37" seat height adjustment
- Foot control standard
- Meets UL 60601-1 and CAN/CSA C22.2 No.601.1-M90 standards
- Table weight: 315 lb
- Shipping weight: 370 lb

Model	Description
223-016	Barrier-Free® Examination Table w/pelvic tilt, drawer heater & receptacle, w/o top
002-0870-XXX ¹	Upholstery top - Soft Touch, Stitched
223-016	Barrier-Free® Examination Table w/pelvic tilt, drawer heater & receptacle, w/o top
002-0871-XXX ¹	Upholstery top - Styled Seamless

¹ - See color code chart on page 30

Midmark Product Guide
Effective June 26, 2017

Accessories	Description
Arm Systems	
9A561001 ^{7,8}	Patient Support Rails, rails only, pair, field installed
9A585001 ⁹	Patient Support Rail Receivers - pair, plant installed
9A314001 ⁴	Adjustable arm system - arms only, field installed
9A314004 ¹⁰	Adjustable arm receiver - plant installed
9A314005	Adjustable arm receiver - field installed
9A354001 ⁵	Assist handles - handles only, 28" wide upholstery only, field installed
9A354002	Assist handles - receiver only, 28" wide upholstery only, plant installed
9A354003	Assist handles - receiver only, 28" wide upholstery only, field installed
Knee Crutches	
9A410001	Standard Knee Crutch only (204/604/22X/62X/230/630) - field installed
9A410004	Standard Knee Crutch Kit with accessory receiver (204/604/22X/62X) - field installed
9A411001	Articulating Knee Crutch only (204/604/22X/62X/230/630) - field installed
9A411004	Articulating Knee Crutch Kit with accessory receiver (204/604/22X/62X) - field installed
9A412003	Mounting Brackets, hardware with accessory receiver (204/604/22X/62X) - field installed
Miscellaneous	
9A104001	Urology drain pan w/hose
9A313001	Rear storage pod, field installed
9A356001	Hand Control
9A363001	PC Work Surface Kit
9A417001	Supply Assistant Storage Unit
9A70000	Stainless steel treatment pan
9A551002	IQvitals® Table Mount, Non-IQhub®
002-0897-02-XXX ^{1,2,6}	Thick footrest pad, stitched (For use with upholstery set 002-0871-XXX)
002-0970-00-XXX ^{1,3,6}	Thick footrest pad, stitched (For use with upholstery set 002-0870-XXX)
002-0971-XXX ¹	Flat stitched upholstery top
002-0971-999	Flat Stitched upholstery top Special Color
002-0870-999	Upholstery top - Soft Touch, stitched Special Color
002-0871-999	Upholstery top - Styled Seamless Special Color
002-0872-XXX ¹	CAL 133 Upholstery Top-Soft Touch Stitched
002-0873-XXX ¹	CAL 133 Upholstery Top-Styled Seamless
002-1573-XXX ¹	CAL 133 Flat stitched upholstery top
002-0872-999	CAL 133 Upholstery Top-Soft Touch Stitched Special Color
002-0873-999	CAL 133 Upholstery Top-Styled Seamless Special Color
002-1573-999	CAL 133 Flat stitched upholstery top Special Color

¹ - See color code chart on page 30

² - For 002-0871-XXX seamless upholstery

³ - For 002-0971-XXX flat stitched upholstery

⁴ - 9A314004 or 9A314005 arm receiver required

⁵ - 9A354003 or 9A354002 assist handle receiver required

⁶ - Additional charge for Special or Premium Upholstery Color for Accessory

⁷ - 9A585001 patient support rail receivers required

⁸ - If also mounting exam light to table, use 253-003 to minimize contact point

⁹ - Cannot order with 9A314004 or 9A354002

¹⁰ - Cannot order with 9A585001



Ritter 222

Barrier-Free® Power Examination Table

- Power height
- Manual back with air spring assist
- 400-pound patient weight capacity
- 18" to 37" seat height adjustment
- Foot control standard
- Meets UL 60601-1 and CAN/CSA C22.2 No.601.1-M90 standards
- Table weight: 295 lb
- Shipping weight: 350 lb

Model	Description
222-016 002-0870-XXX ¹	Barrier-Free® Examination Table w/pelvic tilt, drawer heater & receptacle, w/o top Upholstery top - Soft Touch, Stitched
222-016 002-0871-XXX ¹	Barrier-Free® Examination Table w/pelvic tilt, drawer heater & receptacle, w/o top Upholstery top - Styled Seamless

¹ - See color code chart on page 30

Accessories	Description
Arm Systems	
9A561001 ^{7,8}	Patient Support Rails, rails only, pair, field installed
9A585001 ⁹	Patient Support Rail Receivers - pair, plant installed
9A314001 ⁴	Adjustable arm system - arms only, field installed
9A314004 ¹⁰	Adjustable arm receiver - plant installed
9A314005	Adjustable arm receiver - field installed
9A354001 ⁵	Assist handles - handles only, 28" wide upholstery only, field installed
9A354002	Assist handles - receiver only, 28" wide upholstery only, plant installed
9A354003	Assist handles - receiver only, 28" wide upholstery only, field installed
Knee Crutches	
9A410001	Standard Knee Crutch only (204/604/22X/62X/230/630) - field installed
9A410004	Standard Knee Crutch Kit with accessory receiver (204/604/22X/62X) - field installed
9A411001	Articulating Knee Crutch only (204/604/22X/62X/230/630) - field installed
9A411004	Articulating Knee Crutch Kit with accessory receiver (204/604/22X/62X) - field installed
9A412003	Mounting Brackets, hardware with accessory receiver (204/604/22X/62X) - field installed
Miscellaneous	
9A104001	Urology drain pan w/hose
9A313001	Rear storage pod, field installed
9A363001	PC Work Surface Kit
9A417001	Supply Assistant Storage Unit
9A70000	Stainless steel treatment pan
9A551002	IQvitals® Table Mount, Non-IQhub®
002-0897-02-XXX ^{1,2,6}	Thick footrest pad, stitched (For use with upholstery set 002-0871-XXX)
002-0970-00-XXX ^{1,3,6}	Thick footrest pad, stitched (For use with upholstery set 002-0870-XXX)
002-0971-XXX ¹	Flat stitched upholstery top
002-0872-999	CAL 133 Upholstery Top-Soft Touch Stitched Special Color
002-0873-999	CAL 133 Upholstery Top-Styled Seamless Special Color
002-1573-999	CAL 133 Flat stitched upholstery top Special Color
002-0872-XXX ¹	CAL 133 Upholstery Top-Soft Touch Stitched
002-0873-XXX ¹	CAL 133 Upholstery Top-Styled Seamless
002-1573-XXX ¹	CAL 133 Flat stitched upholstery top
002-0971-999	Flat Stitched upholstery top Special Color
002-0870-999	Upholstery top - Soft Touch, Stitched Special Color
002-0871-999	Upholstery top - Styled Seamless Special Color

¹ - See color code chart on page 30

² - For 002-0871-XXX seamless upholstery

³ - For 002-0971-XXX flat stitched upholstery

⁴ - 9A314004 or 9A314005 arm receiver required

⁵ - 9A354003 or 9A354002 assist handle receiver required

⁶ - Additional charge for Special or Premium Upholstery Color for Accessory

⁷ - 9A585001 patient support rail receivers required

⁸ - If also mounting exam light to table, use 253-003 to minimize contact point

⁹ - Cannot order with 9A314004 or 9A354002

¹⁰ - Cannot order with 9A585001



Ritter 204

Manual Examination Table

- 500-pound patient weight capacity
- Electrical receptacle standard
- Styled Seamless or Soft Touch stitched upholstery
- Pass-thru drawer design, no need to specify RH or LH
- Manual back with air spring assist
- Meets UL 60601-1 and CAN/CSA C22.2 No.601.1-M90 standards
- Please specify upholstery/drawer panel color
- Bulk storage models have Pebble Grey doors
- Table weight: 250 lb
- Shipping weight: 300 lb

Model	Description
Pass-Thru Drawer Models	
204-001 002-0870-XXX ¹	Manual Examination Table w/o top Upholstery top - Soft Touch, Stitched
204-001 002-0871-XXX ¹	Manual Examination Table w/o top Upholstery top - Styled, Seamless
204-002 002-0870-XXX ¹	Manual Examination Table w/drawer heater, pelvic tilt w/o top Upholstery top - Soft Touch, Stitched
204-002 002-0871-XXX ¹	Manual Examination Table w/drawer heater, pelvic tilt w/o top Upholstery top - Styled, Seamless
Bulk Storage Models	
204-003 002-0870-XXX ¹	Manual Examination Table, Bulk Storage, w/o top Upholstery top - Soft Touch, Stitched
204-003 002-0871-XXX ¹	Manual Examination Table, Bulk Storage, w/o top Upholstery top - Styled, Seamless
204-004 002-0870-XXX ¹	Manual Examination Table, Bulk Storage w/drawer heater & pelvic tilt w/o top Upholstery top - Soft Touch, Stitched
204-004 002-0871-XXX ¹	Manual Examination Table, Bulk Storage w/drawer heater & pelvic tilt w/o top Upholstery top - Styled, Seamless

¹ - See color code chart on page 30

Midmark Product Guide

Effective June 26, 2017

Accessories	Description
Arm/Handle Systems	
9A314001 ²	Adjustable arm system (arms only)
9A354001 ²	Assist handles (handles only), 28" wide upholstery only, field installed
9A354002	Accessory receiver only, 28" wide upholstery only, plant installed
9A354003	Accessory receiver only, 28" wide upholstery only, field installed
Knee Crutches	
9A410001	Standard Knee Crutch only (204/604/22X/62X/230/630) - field installed
9A410004	Standard Knee Crutch Kit with accessory receiver (204/604/22X/62X) - field installed
9A411001	Articulating Knee Crutch only (204/604/22X/62X/230/630) - field installed
9A411004	Articulating Knee Crutch Kit with accessory receiver (204/604/22X/62X) - field installed
9A412003	Mounting Brackets, hardware with accessory receiver (204/604/22X/62X) - field installed
Miscellaneous	
9A104001	Urology drain pan w/hose
9A363001	PC Work Surface Kit
9A372001	Pass-through drawer stop (set of 3)
9A70000	Stainless steel treatment pan
9A551002	IQvitals® Table Mount, Non-IQhub®
002-0965-00-216	Pebble Grey drawer inserts (set of 8)
002-0897-02-XXX ^{1,3}	Thick footrest pad, stitched (For use with upholstery set 002-0871-XXX)
002-0971-XXX ¹	Flat Stitched Upholstery Top
002-0971-999	Flat Stitched upholstery top Special Color
002-0870-999	Upholstery top - Soft Touch, Stitched Special Color
002-0871-999	Upholstery top - Styled Seamless Special Color
002-0970-XXX ^{1,3}	Thick footrest pad, stitched (For use with upholstery set 002-0870-XXX)
002-0872-XXX ¹	CAL 133 Upholstery Top-Soft Touch Stitched
002-0873-XXX ¹	CAL 133 Upholstery Top-Styled Seamless
002-1573-XXX ¹	CAL 133 Flat stitched upholstery top
002-0872-999	CAL 133 Upholstery Top-Soft Touch Stitched Special Color
002-0873-999	CAL 133 Upholstery Top-Styled Seamless Special Color
002-1573-999	CAL 133 Flat stitched upholstery top Special Color

¹ - See color code chart on page 30

² - 9A354002 or 9A354003 accessory receiver required

³ - Additional charge for Special or Premium Upholstery Color for Accessory

Ritter 203

Treatment Table



- 400-pound patient weight capacity
- 27" W x 74" D x 31" H
- Please specify upholstery /panel color
- Panels match upholstery color
- Table weight: 135 lb
- Shipping weight: 170 lb

Model	Description
203-001-XXX ¹	Treatment Table Standard
203-002-XXX ¹	Treatment Table w/drawers
203-003*-XXX ¹	Treatment Table w/stirrups
203-004*-XXX ¹	Treatment Table w/drawers and stirrups

* Stirrups factory installed only
¹ - See color code chart on page 30

Miscellaneous	Description
203-001-999	Treatment Table Standard - Special Color
203-002-999	Treatment Table Standard - Special Color
203-003-999	Treatment Table Standard - Special Color
203-004-999	Treatment Table Standard - Special Color
203-005-XXX ¹	CAL 133 Treatment Table Standard
203-006-XXX ¹	CAL 133 Treatment Table w/drawers
203-007-XXX ¹	CAL 133 Treatment Table w/stirrups
203-008-XXX ¹	CAL 133 Treatment Table w/drawers and stirrups
203-005-999	CAL 133 Treatment Table Standard - Special Color
203-006-999	CAL 133 Treatment Table w/drawers - Special Color
203-007-999	CAL 133 Treatment Table w/stirrups - Special Color
203-008-999	CAL 133 Treatment Table w/drawers and stirrups - Special Color

Options	Description
AD 1468	Table less pillow (-001: -002 models)
AD 1469	Table less pillow (-003: -004 models)
AD 1470	Special height: 26", 27", 28", & 29" only

Ritter 95

Treatment Table



- 350-pound patient weight capacity
- 27" W x 72" D x 30" H
- Please specify upholstery color
- Frame is a dark oak finish only
- Paper roll holder included
- Table weight: 80 lb
- Shipping weight: 100 lb


Model	Description
95-001-XXX ¹	Treatment Table w/o shelf
95-001-999	Treatment Table w/o shelf Special Color

Accessories	Description
9A364001-XXX ¹	Shelf (field installed only)
9A364001-999	Shelf Special Color

¹ - See color code chart on page 30

Options	Description
95-002-XXX ¹	CAL 133 Treatment Table w/o shelf
95-002-999	CAL 133 Treatment Table w/o shelf - Special Color

Ritter Color-Specific Product Number Chart



Example:

204-001 *Table*

002-0870-232 *Top*

204 Examination Table - Shadow, Soft Touch, Stitched Upholstery Top & Panels

216	Pebble Grey
229	Perfect Plum
230	Moss
231	Navy
232	Shadow
233	Dusty Blue
234	Clay
312	Black



Midmark 425

Air Lift Physician Stool (Hand Operated)

- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Model

425-001
002-0273-XXX^{1,2}

Description

Basic stool w/o seat cushion
Seat cushion

Accessories

002-0274-XXX^{1,3}
002-0275-XXX^{1,3}
002-1546-00

Back kit
Arm kit
Lab height w/black foot ring (425 only)

Miscellaneous

002-0944-XXX¹
002-0945-XXX¹
002-0946-XXX¹

CAL 133 Seat Cushion-Premium Upholstery
CAL 133 Back Kit-Premium Upholstery
CAL 133 Arm Kit-Premium Upholstery

¹ - See color code chart on page 36

² - Additional charge for Special or Ultra Premium Upholstery Color

³ - Additional charge for Special or Ultra-Premium Upholstery Color for Accessory



Midmark 427

Air Lift Physician Stool (Foot Operated)

- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Model

427-001
002-0273-XXX^{1,2}

Description

Basic stool w/o seat cushion
Seat cushion

Accessories

002-0274-XXX^{1,3}
002-0275-XXX^{1,3}

Back kit
Arm kit

Miscellaneous

002-0944-XXX¹
002-0945-XXX¹
002-0946-XXX¹

CAL 133 Seat Cushion-Premium Upholstery
CAL 133 Back Kit-Premium Upholstery
CAL 133 Arm Kit-Premium Upholstery

¹ - See color code chart on page 36

² - Additional charge for Special or Ultra Premium Upholstery Color

³ - Additional charge for Special or Ultra-Premium Upholstery Color for Accessory



Ritter 279

Air Lift Surgeon Stool (Hand Operated)

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 30 lb
- Shipping weight: 35 lb

Model

279-001-XXX^{1,2}

Description

Basic stool w/hand release - pneumatically adjustable



Ritter 278

Air Lift Surgeon Stool (Foot Operated)

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 30 lb
- Shipping weight: 35 lb

Model

278-001-XXX^{1,2}

Description

Basic stool w/foot release - pneumatically adjustable



hand release

Ritter 277

Air Lift Stool with Backrest

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Model

277-001-XXX^{1,2}

277-201-XXX^{1,2}

Description

Basic stool w/back - pneumatically adjustable, hand release

Basic stool w/back - pneumatically adjustable, foot release



foot release



hand release

Ritter 276

Air Lift Stool

- Chrome five leg base
- Weight capacity: 350 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Model

276-001-XXX^{1,2}

276-201-XXX^{1,2}

Description

Basic stool - pneumatically adjustable, hand release

Basic stool - pneumatically adjustable, foot release



foot release

¹ - See color code chart on page 36

² - Additional charge for Special or Ultra-Premium Upholstery Color



Ritter 275

Adjustable Physician Stool with Backrest

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Model
275-001-XXX^{1,2}

Description
Basic stool w/back - manually adjustable



Ritter 274

Adjustable Physician Stool

- Chrome five leg base
- Weight capacity: 350 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Model
274-001-XXX^{1,2}

Description
Basic stool - manually adjustable



Ritter 273

Air Lift Stool with Backrest

- Black composite five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Model
273-001-XXX^{1,2}

Description
Basic stool w/back - pneumatically adjustable



Ritter 272

Air Lift Stool

- Black composite five leg base
- Weight capacity: 300 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Model
272-001-XXX^{1,2}

Description
Basic stool - pneumatically adjustable

¹ - See color code chart on page 36

² - Additional charge for Special or Ultra-Premium Upholstery Color



Model
271-001-XXX^{1,2}

Ritter 271

Adjustable Stool with Backrest

- Black composite five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description
Basic stool w/back - manually adjustable



Model
270-001-XXX^{1,2}

Ritter 270

Adjustable Stool

- Black composite five leg base
- Weight capacity: 300 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Description
Basic Stool - manually adjustable

Model

Miscellaneous

27X-X11-XXX¹

9A366001-XXX^{1,3}

9A366002-XXX¹

27X-002-XXX^{1,2}

27X-003-XXX^{1,2}

27X-004-XXX^{1,2}

27X-005-XXX^{1,2}

9A367001

9A368001

Description

CAL 133 Seat Upholstery

Back Rest Kit (Weight capacity, 225 lb)

CAL 133 Back Rest Kit (Weight capacity, 225 lb)

Set of Soft Rubber Casters

Set of Glides (not available on foot release models)

Set of Locking Casters

Auto Locking Casters

Foot ring, field installed (not available on foot release models)

6" Lab Height Ext. Kit (not available on foot release models)

¹ - See color code chart on page 36

² - Additional charge for Special or Ultra-Premium Upholstery Color

³ - Additional charge for Special or Ultra-Premium Upholstery Color for Accessory



Midmark 680

Upholstered Side Chair

- 400-pound patient weight capacity
- Wall saver, kickback rear leg design
- Leveling feet
- Chair weight: 20 lb
- Shipping weight: 30 lb

Model	Description
680-001-XXX ¹	Basic chair (without arms)
680-001-XXX ¹	Basic chair (without arms) Ultra-Premium Upholstery
680-002-XXX ¹	Basic chair (with arms)
680-002-XXX ¹	Basic chair (with arms) Ultra-Premium Upholstery

Miscellaneous

680-001-999	Basic chair (without arms) Special Color
680-002-999	Basic chair (with arms) Special Color
680-003-XXX ¹	CAL 133 Basic Chair (without arms)
680-004-XXX ¹	CAL 133 Basic Chair (with arms)
680-003-999	CAL 133 Basic Chair (without arms) Special Color
680-004-999	CAL 133 Basic Chair (with arms) Special Color

305¹ - See color code chart on page 36



Ritter 281

Blood Drawing Chair

- 400-pound patient weight capacity
- Optional storage drawer
- Chair weight: 60 lb
- Shipping weight: 75 lb



Model	Description
281-001-XXX ¹	Basic chair (without storage drawer)
281-001-XXX ¹	Basic chair (without storage drawer) Ultra-Premium Upholstery
281-002-XXX ¹	Basic chair (with storage drawer)
281-002-XXX ¹	Basic chair (with storage drawer) Ultra-Premium Upholstery


Miscellaneous

281-001-999	Basic chair (without storage drawer) Special Color
281-002-999	Basic chair (with storage drawer) Special Color
281-003-XXX ¹	CAL 133 Basic Chair (without storage drawer)
281-004-XXX ¹	CAL 133 Basic Chair (with storage drawer)
281-003-999	CAL 133 Basic Chair (without storage drawer) Special Color
281-004-999	CAL 133 Basic Chair (with storage drawer) Special Color

¹ - See color code chart on page 36

Seating Color-Specific Product Number Chart

XXX - XXX - XXX



Example: 276-001-230
276 Air Lift Stool, Moss Upholstery

Midmark
Solid Premium Colors

235	Blueberry
236	Silver Sage
237	Tea Green
238	Fossil Grey
239	Cashmere Blue
240	Terra Cotta
312	Black

Midmark
Pattern Premium Colors

247	Belagio
248	Siena
249	Tuscany
250	Firenze
251	Capri
252	Treviso

Midmark
Ultra-Premium Colors¹

228	Ultra-Raven Wing
241	Ultra-Pacific Blue
242	Ultra-Spruce
243	Ultra-Olivine
244	Ultra-Arctic
245	Ultra-Sky
246	Ultra-Canyon

¹ See product price sheet
for Ultra-Premium color
availability and price.

Ritter Solid Colors

216	Pebble Grey
229	Perfect Plum
230	Moss
231	Navy
232	Shadow
233	Dusty Blue
234	Clay
312	Black



Mobile Procedure Light



Single Procedure Light



Dual Procedure Light



Wall Mount

Ritter 255

LED Procedure Light

- 4,400K color temperature
- Color Rendering Index (CRI) greater than 95
- 7,500 fc at 1 meter (39 inches) focal length
- Effortless positioning
- IEC/ES 60601-1, CAN/CSA C22.2 60601-1, 60601-2-41
- Weight: 65 lb (118 lb for dual light)
- Shipping weight: 120 lb (190 lb for dual light)

Model

Description

255-001	Single mount, 8' ceiling
255-002	Single mount, 9' ceiling
255-003	Dual mount, 8' ceiling
255-004	Dual mount, 9' ceiling
255-005	Mobile Procedure Light
255-006	Wall Mount

Accessories

002-1865-00 ^{2,3}	10' ceiling mount adapter kit
9A584001 ¹	Conversion kit w/cover (single to double light)

¹ - To convert an existing single light to a dual light, this conversion kit must be ordered in addition to a second 8' or 9' light

² - Converts an 8' or 9' light to a 10' ceiling height

³ - One required per light arm

The following model numbers are recommended based on the following ceiling heights (measured from the bottom of the ceiling surface to the floor):

Less than 8'	255-005	Mobile Procedure Light
8' to 9'	255-001, 255-003	8' Single or Dual Procedure Light
9' to 10'	255-002, 255-004	9' Single or Dual Procedure Light
Greater than 10'	255-001, 255-003, 255-005	8' Single or Dual Light with 10' ceiling mount adapter kit(s) or Mobile Procedure Light

Ritter 253

LED Exam Light



- LED technology for up to 50,000 hours of useful life
- Designed to consume 45-60% less energy than other lights
- 4,500K color temperature
- 1,700 fc intensity (18,298 lux) at 16-inch focal length
- Sleek, ergonomic and easy-to-clean design
- Weight of light unit: 6 lb (253-002 light kit including bracket)
- Shipping weight: 8 lb (253-002 light kit including bracket)
- UL6060-1 and CAN/CSA C22.2 No. 601.1-M90



Exam Table Kit

Mobile Caster Stand Kit

Model	Description
253-001	LED Exam Light - light only - excludes mounting hardware
253-002 ³	LED Exam Light - exam table kit - models 204/22X*/604/62X* - includes light, splines and table mounting hardware, field installed
253-003 ²	LED Exam Light - exam table kit - model 625 for 32" upholstery and Patient Support Rails - includes light, splines and table mounting hardware, field installed
253-005	LED Exam Light - procedure chair kit - models 230/641/646/647 - includes light, base rail and chair mounting hardware, field installed
253-006	LED Exam Light - mobile caster stand kit - includes light & caster stand
253-007	LED Exam Light - wall mount kit - includes light, wall mounting bracket and hardware, field installed

Miscellaneous

9A454001	LED exam table mounting hardware (only) - models 204/22X*/604/62X* - field installed
9A455001	LED exam table mounting hardware (only) - model 625 w/ 32" upholstery - field installed
9A457001	LED procedure chair mounting hardware (only) - models 230/641/646/647 - field installed
9A562001	LED 630 procedure table spline mounting hardware (only), 28" wide upholstery, field installed
9A563001	LED 630 procedure table spline mounting hardware (only), 32" wide upholstery, field installed
9A564001 ¹	LED 630 procedure table rail mount hardware (only), field installed
9A459001	Disposable Sterility Sleeve; 5 boxes of 25 per carton
9A458001	LED light mobile caster stand (only)
9A474001	LED wall mount bracket (only)

* Mounts to tables with a 2006 or newer manufacture date

¹ - Requires (1) accessory rails 9A557001 or 9A557002

² - Use 253-003 with Patient Support Rails

³ - Slight interference with Patient Support Rails, use 253-003 to minimize



Ritter 250

LED Exam Light

- LED technology for up to 25,000 hours of useful life
- Designed to consume 45-60% less energy than other lights
- 4,000K color temperature
- 1,800 fc intensity (19,375 lux)
- UL6060-1 and CAN/CSA C22.2 No. 601.1-M90
- Exam light weight: 20 lb
- Shipping weight: 25 lb



Caster Base Accessory

Model	Description
250-002	LED Exam Light
9A274001	Caster base accessory

Midmark Diagnostic Medical Products (ECG, Holter, Spirometer, Stress, and Vital Signs) Limited Warranty

SCOPE OF WARRANTY

Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will repair or replace IQ brand devices manufactured by Midmark (except for those not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of applicable devices. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which, upon examination by Midmark, prove to be defective. This warranty extends only to the original retail purchaser of a product, and is not transferable or assignable.

Midmark warrants to the original retail purchaser that it will repair or replace software contained within the IQ brand products manufactured by Midmark (except for those not warranted under "Exclusions") if: (1) the media on which the software is furnished exhibits defects in material or workmanship under normal use; or (2) the software does not substantially conform to its published specifications. Midmark does not warrant that the software: (1) is error free; (2) can be used without problems or interruptions; or (3) is free from vulnerability to intrusion or attack by viruses or other methods. Midmark also agrees to provide the original retail purchaser with software updates, if any, that become generally available to the public within the warranty period stated, for the specific product owned by the purchaser.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original retail purchaser, shall be as follows:

- (a) One (1) year for software and diagnostic devices (unless otherwise noted below and except items listed under "Exclusions"),
- (b) Two (2) years for vital signs devices,
- (c) One (1) year for reusable blood pressure cuffs and SpO₂ sensors,
- (d) Ninety (90) days for all other accessories (sold in a product kit or purchased separately),

OBTAINING WARRANTY SERVICE

Warranty service must be obtained through Midmark. Contact Midmark Technical Service for warranty service inquiries or issues by phone at 1-800-624-8950, option 2, or via email at techsupport@midmark.com.

It is the original retail purchaser's obligation to arrange for delivery of a product to Midmark for warranty service or exchange. It is also the original retail purchaser's obligation to comply with the warranty service instructions provided by Midmark.

EXCLUSIONS

This warranty does not cover and Midmark shall not be liable for the following:

- (1) any product that has been repaired, modified, or tampered with by any person or entity other than Midmark;
- (2) any product for which any label or tag attached to the product has been removed or tampered with;
- (3) defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering, or failure to seek and obtain repair or replacement in a timely manner;
- (4) products which are not installed, used, and properly cleaned and maintained as stated in the Midmark Installation, User Guide and/or Operation Manual for the applicable product;
- (5) products considered to be of a consumable nature;
- (6) accessories or parts covered by the respective original manufacturers' warranty, including but not limited to treadmills, computers, printers, scales, carts, mounting options, and carrying cases, after the initial 90 days;
- (7) charges by anyone for adjustments, repairs, replacement parts, installation, or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (8) costs and expenses of routine maintenance and cleaning;
- (9) representations and warranties made by any person or entity other than Midmark; and
- (10) items provided through sales promotions.

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER

MIDMARK'S ONLY OBLIGATION UNDER THIS WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF DATA, LOSS OF USE, DOWNTIME, COVER, AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS, AND BENEFITS.

NO AUTHORIZATION

No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

WARRANTY DISCLAIMER

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

THIS WARRANTY WILL BE VOID IF ANY COVERED PRODUCT IS REPAIRED, MODIFIED, OR TAMPERED WITH BY ANY PERSON OR ENTITY OTHER THAN MIDMARK, OR IF ANY LABEL OR TAG ATTACHED TO A COVERED PRODUCT IS REMOVED OR TEMPERED WITH.

STATUTE OF LIMITATIONS

No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period. In the event multiple warranty periods exist with respect to a product, the ninety (90) day period provided for herein shall begin to run from expiration of the warranty period for the component to which the claim relates.

Midmark Diagnostic Medical Products Support Information

Customer Support Services

For technical assistance with a Midmark diagnostic medical product, please contact Technical Service at (800) 624-8950, option 2, Monday through Friday, 6:00 am to 4:00 pm, Pacific time. For troubleshooting tips and product use information refer to the operation manual. Additional help topics, self-guided user tools and documentation can be accessed at any time through our Support section and knowledge base at midmark.com.

Return Material Authorization

To return any product for repair, an RMA (Return Material Authorization) number must be obtained from a Midmark Technical Service Specialist. The RMA number should be referenced on the package(s) containing the items to be returned.

Shipping

Before shipping any product to Midmark, be certain that an RMA (Return Material Authorization) number has been issued and that all guidelines regarding this authorization are followed. We highly recommend that you follow all guidelines for the shipment of medical products set forth by the shipping company you choose to use. If a question should arise regarding the appropriate method of shipment, please feel free to ask when calling for your RMA number. It is the responsibility of the customer when shipping a product to ensure that all packages and their contents get to Midmark safely.

Midmark will not assume responsibility for damage incurred by improper packaging, shipment or product use. Such actions will void all applicable warranties.

Midmark IQ Devices Part Number Ordering Information

Note: Midmark recommends that all Midmark IQ devices drop ship directly to the customer

TO ORDER:

Midmark IQ devices require either an EMR that has integrated Midmark software or Midmark IQmanager® Software to operate. Devices include, but are not limited to, the part numbers shown below:

4-000-0025
4-000-0026
4-000-0062
4-000-0110
4-000-0113
4-000-0116

If the customer **currently has** an Electronic Medical Record (EMR) software package that has integrated Midmark Software, no additional software is required.

Place an order for the appropriate device(s) and please provide:

- EMR name and version
- Customer facility & contact name, address and phone
- Customer computer environment (Thin Client, Fat Client)

If the customer **does not have** an EMR, or using Allscripts TouchWorks™ EHR, Epic® EMR or GE Centricity®: please contact your local Midmark Sales Representative to order the appropriate Midmark software solution. Additional charges or license fees will apply.

For ordering questions regarding Midmark IQ devices, please contact your local Midmark Sales Representative or call 800-624-8950, option 2.

Midmark IQvitals® Zone™ Vital Signs Wireless Monitor

with Touchscreen, Serenity™ BP and Exergen TemporalScanner™

- Includes Spot, Averaging, Continuous and Interval BP Modes



Shown with IQvitals® Zone™
Mobile Cart Option

Mobile cart not included in pricing

- 1 - IQvitals® Zone™ Wireless Monitor
- 1 - Reusable Dual Tube BP Cuff – Adult (26-35cm)
- 1 - Reusable Dual Tube BP Cuff – Large Adult (32-42cm)
- 1 - Dual Tube BP Hose – 6.5'
- 1 - Exergen TemporalScanner™ Thermometer
- 1 - Exergen Disposable Caps (10 ea)
- 1 - IQvitals® Zone™ Wireless USB Dongle
- 1 - IQvitals® Zone™ USB Cable – 10'
- 1 - IQvitals® Zone™ Laptop Tray – 15"
- 1 - Front Storage Bin Kit
- 1 - Device to Temperature Module Cable
- 1 - Three-Arm Knobs (Set of 2)
- 1 - AC Power Supply
- 1 - Setup Guide
- 1 - Quick Reference Guide
- 1 - Operation Manual CD

Model	Description
4-000-0724	IQvitals® Zone™, Wireless Monitor, with Serenity™ BP, Exergen TemporalScanner™ and Masimo SpO ₂
4-000-0722	IQvitals® Zone™, Wireless Monitor, with Serenity™ BP, Exergen TemporalScanner™ and Nellcor™ SpO ₂
4-000-0720	IQvitals® Zone™, Wireless Monitor, with Serenity™ BP and Exergen TemporalScanner™

Midmark IQvitals® Zone™ Non-Wireless Vital Signs Device

with Touchscreen, Serenity™ BP and Exergen TemporalScanner™

- Includes Spot and Averaging BP Modes



Shown with IQvitals® Zone™
Mobile Cart Option

- 1 - IQvitals® Zone™ Vital Signs Device
- 1 - Reusable Dual Tube BP Cuff – Adult (26-35cm)
- 1 - Reusable Dual Tube BP Cuff – Large Adult (32-42cm)
- 1 - Dual Tube BP Hose – 6.5'
- 1 - Exergen TemporalScanner™ Thermometer
- 1 - Exergen Disposable Caps (10 ea)
- 1 - IQvitals® Zone™ USB Cable – 10'
- 1 - IQvitals® Zone™ Laptop Tray – 15"
- 1 - Front Storage Bin Kit
- 1 - Device to Temperature Module Cable
- 1 - Three-Arm Knobs (Set of 2)
- 1 - AC Power Supply
- 1 - Setup Guide
- 1 - Quick Reference Guide
- 1 - Operation Manual CD

Mobile cart not included in pricing

Model	Description
4-000-0624	IQvitals® Zone™, non-Wireless, with Serenity™ BP, Exergen TemporalScanner™ and Masimo SpO ₂
4-000-0622	IQvitals® Zone™, non-Wireless, with Serenity™ BP, Exergen TemporalScanner™ and Nellcor™ SpO ₂
4-000-0620	IQvitals® Zone™, non-Wireless, with Serenity™ BP and Exergen TemporalScanner™

IQvitals® Zone™ Mounting Options

Model	Description
2-200-0072	IQvitals® Zone™ Mobile Cart
2-200-0082	IQvitals® Zone™ Wall Mount Articulating Arm

IQvitals® Zone™ Accessories and Supplies

Model	Description
1-100-1603	Fairbanks® TeleWeigh® Digital Scale
3-200-0030	Blood Pressure Cuff, Dual Tube, Reusable, Child (13-20cm)
3-200-0031	Blood Pressure Cuff, Dual Tube, Reusable, Small Adult (18-26cm)
3-200-0032	Blood Pressure Cuff, Dual Tube, Reusable, Adult (26-35cm)
3-200-0033	Blood Pressure Cuff, Dual Tube, Reusable, Adult Long (29-38cm)
3-200-0034	Blood Pressure Cuff, Dual Tube, Reusable, Large Adult (32-42cm)
3-200-0035	Blood Pressure Cuff, Dual Tube, Reusable, Large Adult Long (35-44cm)
3-200-0036	Blood Pressure Cuff, Dual Tube, Reusable, Thigh (42-50cm)
2-200-0020	Blood Pressure Cuff, Dual Tube, Disposable, Child (13-20cm), 5/pk
2-200-0021	Blood Pressure Cuff, Dual Tube, Disposable, Small Adult (18-26cm), 5/pk
2-200-0022	Blood Pressure Cuff, Dual Tube, Disposable, Adult (26-35cm), 5/pk
2-200-0023	Blood Pressure Cuff, Dual Tube, Disposable, Adult Long (29-38cm) , 5/pk
2-200-0024	Blood Pressure Cuff, Dual Tube, Disposable, Large Adult (32-42cm), 5/pk
2-200-0025	Blood Pressure Cuff, Dual Tube, Disposable, Large Adult Long (35-44cm), 5/pk
2-200-0026	Blood Pressure Cuff, Dual Tube, Disposable, Thigh (42-50cm), 5/pk
2-200-0030	SpO ₂ Sensor, Reusable, Masimo, Pediatric (M-LNCS DCI-P)
2-200-0032	SpO ₂ Sensor, Reusable, Masimo, Adult (M-LNCS DC-I)
2-200-0031	SpO ₂ Patient Cable, Masimo (RC-4) - 4'
3-200-0007	SpO ₂ Sensor Kit, Reusable, Nellcor, Pediatric (D-YS and D-YSPD)
2-200-0233	SpO ₂ Sensor, Reusable, Nellcor, Adult (DS100A)
2-200-0234	SpO ₂ Interface Cable, Nellcor (DOC4) - 4'
2-200-0115	Disposable Caps, Exergen Thermometer, 1000pc/box
2-200-0116	Disposable Caps, Exergen Thermometer, 100pc/box
2-200-0074	Storage Bin Kit, Dual, Cart Mounted
2-200-0084	Storage Bin Kit, Wall Mounted
3-200-0009	Tray, IQvitals® Zone™ - 12"
3-200-0012	USB Dongle, IQvitals® Zone™ Wireless

Connectivity Supplies for Midmark 625

Model	Description
9A478001	Midmark IQscale® serial cable for IQvitals®, 6 foot length, coiled
9A478002	Midmark IQscale® serial cable for IQvitals®, 15 foot length, straight
9A478003	Midmark IQscale® serial cable for IQvitals®, 30 foot length, straight
9A478004	Midmark IQscale® serial cable for IQvitals®, 50 foot length, straight
002-10063-00	Connectivity Wiring Kit (For 625-004)

Midmark IQvitals®



Vital Signs Device with Touchscreen - Blood Pressure and Temperature

- 1 - IQvitals® Device w/Blood Pressure and Temperature
- 1 - AC Power Supply
- 1 - Reusable Blood Pressure Cuff - Adult
- 1 - Reusable Blood Pressure Cuff - Large Adult
- 1 - Blood Pressure Hose - 6.5'
- 1 - Oral/Axillary Temperature Probe
- 1 - Temperature Probe Covers (box of 20)
- 1 - Serial Cable - 10'
- 1 - USB Cable - 10'
- 1 - Operation Manual CD
- 1 - Quick Reference Guide

Model

4-000-0500

Description

IQvitals® with Blood Pressure and Alaris® Turbo Temp® Thermometer

Vital Signs Device with Touchscreen - Blood Pressure, Temperature and SpO₂

- 1 - IQvitals® Device w/Blood Pressure, Temperature and SpO₂
- 1 - AC Power Supply
- 1 - Reusable Blood Pressure Cuff - Adult
- 1 - Reusable Blood Pressure Cuff - Large Adult
- 1 - Blood Pressure Hose - 6.5'
- 1 - Reusable SpO₂ Adult Finger Sensor
- 1 - SpO₂ Extender Cable - 4'
- 1 - Oral/Axillary Temperature Probe
- 1 - Temperature Probe Covers (box of 20)
- 1 - Serial Cable - 10'
- 1 - USB Cable - 10'
- 1 - Operation Manual CD
- 1 - Quick Reference Guide

Model

4-000-0510

Description

IQvitals® with Blood Pressure, Alaris® Turbo Temp® Thermometer and SpO₂

IQvitals® Accessories and Supplies

Model	Description
1-100-1603	Fairbanks® TeleWeigh® Digital Scale
1-100-1605	IQvitals® Thermal Printer
3-004-2004	IQvitals® Mobile Cart Printer Bracket
3-004-2010	IQvitals® Mobile Cart Scale Mount
3-009-0010	AC Power Supply
3-009-0012	Power Supply Adapter Kit
3-009-0014	IQvitals® & IQvitals® PC Replacement lithium ion battery
3-009-0016	USB Cable - 10'
3-009-0018	Serial Cable - 10'
3-009-0022	Blood Pressure Hose - 6.5'
3-009-0100	Blood Pressure Hose - 10'
3-009-0021	Reusable SpO ₂ Small Finger Sensor, Envitec®
3-009-0020	Reusable SpO ₂ Adult Finger Sensor, Envitec®
3-009-0026	SpO ₂ Extender Cable, 4', Envitec®
3-009-0024	Alaris® Turbo Temp® Oral/Axillary Temperature Probe
3-009-0058	Turbo Temp® Probe Covers, Carton of 10 boxes (20 covers per box)
3-009-0068	Reusable Blood Pressure Cuff, Infant (8-14 cm)
3-009-0070	Reusable Blood Pressure Cuff, Child (13-20 cm)
3-009-0062	Reusable Blood Pressure Cuff, Small Adult (18-26 cm)
3-009-0064	Reusable Blood Pressure Cuff, Adult, Midmark (26-35 cm)
3-009-0066	Reusable Blood Pressure Cuff, Large Adult, Midmark (32-42 cm)
3-009-0072	Reusable Blood Pressure Cuff, Adult Long (29-38 cm)
3-009-0074	Reusable Blood Pressure Cuff, Large Adult Long (35-44 cm)
3-009-0076	Reusable Blood Pressure Cuff, Thigh (42-50 cm)
3-004-2000	IQvitals® Mobile Cart
3-009-0001	IQvitals® Countertop Mount
3-009-0003	IQvitals® Wall Mount with Basket
3-004-2008	IQvitals® Equipment Pole Mount

Connectivity Supplies for Midmark 625

Model	Description
9A551001	IQvitals® Table Mount, IQhub®
9A551002	IQvitals® Table Mount, Non-IQhub®
9A478001	Midmark IQscale® serial cable for IQvitals®, 6 foot length, coiled
9A478002	Midmark IQscale® serial cable for IQvitals®, 15 foot length, straight
9A478003	Midmark IQscale® serial cable for IQvitals®, 30 foot length, straight
9A478004	Midmark IQscale® serial cable for IQvitals®, 50 foot length, straight
002-10063-00	Connectivity Wiring Kit (For 625-004)



Midmark IQspiro® Digital Spirometer

- 1 - Midmark IQspiro® Digital Spirometer
- 10 - Disposable Spirometer Mouthpieces
- 1 - 10/pack Disposable Nose Clips
- 1 - Mouse Pad
- 1 - Operation Manual CD
- 1 - Quick Reference Guide
- 1 - Carrying Case

Model	Description
4-000-0026	Midmark IQspiro® - USB version with syringe
4-000-0025	Midmark IQspiro® - USB version

Software

Please contact your local Midmark Sales Representative to order the appropriate Midmark software solution, including the new IQmanager® version 10.0.0 for non-EMR users.

IQspiro® Accessories and Supplies

Model	Description
1-100-0007	IQcal® 3-Liter Calibration Syringe
2-100-1205	Disposable Spirometer Mouthpieces (box of 25)
2-100-1206	Disposable Spirometer Mouthpieces (case - 4 boxes)
2-100-1207	Box of 500 Mouthpieces, individually bagged, bulk
2-244-0001	Disposable Nose Clips (box of 25)



Midmark IQecg® Digital ECG with Lead Management (Interpretive)

- 1 - Midmark IQecg® Digital ECG Acquisition Module
- 1 - 10-Lead Patient Cable, with Lead Management
- 1 - 10/pack Clear ECG Clips
- 1 - 100/pack Disposable EZ-Trode ECG Electrodes
- 1 - Mouse Pad
- 1 - Operation Manual CD
- 1 - Quick Reference Guide
- 1 - Carrying Case

Model	Description
4-000-0062	Midmark IQecg® - with Lead Management

Software

Please contact your local Midmark Sales Representative to order the appropriate Midmark software solution, including the new IQmanager® version 10.0.0 for non-EMR users.

IQecg® Supplies

Model	Description
3-047-0001	Universal ECG Clips (3mm, 4mm & Snap) 10/pack
3-047-0005	Clear ECG Clips (4mm & Snap) 10/pack
3-100-0203	Patient Cable, with Lead Management, IQecg®
2-100-0205	Disposable EZ-Trode ECG Electrodes (box of 500)
2-100-0206	Disposable EZ-Trode ECG Electrodes (case - 10 boxes)

Midmark IQholter® Digital Holter



Holter Kit with Digital Recorder

- 1 - Midmark IQholter® Recorder
 - Includes: lead wire set, secure digital card
- 1 - Reusable IQholter® Recorder Pouch and Belt
- 1 - Secure Digital Card Reader
- 1 - Software Security Key, USB Version
- 4 - Holter Prep Kits
- 1 - Mouse Pad
- 1 - Operation Manual CD
- 1 - Quick Reference Guide

Model	Description
4-000-0110	Midmark IQholter® with one (1) Recorder
4-000-0113	Midmark IQholter® EX with one (1) Recorder
4-000-0116	Midmark IQholter® EP with Pacemaker detection and one (1) Recorder

Digital Holter Recorder

- 1 - Midmark IQholter® Recorder
- 1 - Lead wire set
- 1 - Secure Digital Card
- 1 - Operation Manual CD
- 4 - Holter Prep Kits
- 1 - Reusable IQholter® Recorder Pouch and Belt

Model	Description
4-000-0119	Midmark IQholter® Recorder
1-370-0020	Midmark IQholter® Recorder Upgrade Kit (24, 48, 72 Hours) One Upgrade Kit Required per Recorder

Software

Please contact your local Midmark Sales Representative to order the appropriate Midmark software solution, including the new IQmanager® version 10.0.0 for non-EMR users.

Secure Memory Card

Model	Description
1-370-0010	Midmark IQholter® Secure Digital Memory Card
1-370-0004	USB Secure Digital Card Reader

IQholter®, EX, EP Digital Holter Supplies

Model	Description
2-100-0080	Prep Kit - Midmark IQholter®, 5-lead (each) Each kit contains: 1 diary, 2 AA batteries, 1 razor 1 Abrading Pad, 2 alcohol wipes, 5 electrodes
2-100-0081	Prep Kit - Midmark IQholter®, 5-lead (case) Each kit contains: 1 diary, 2 AA batteries, 1 razor 1 Abrading Pad, 2 alcohol wipes, 5 electrodes (case - 30 kits)
2-100-0090	Prep Kit - Midmark IQholter®, 7-lead (each) Each kit contains: 1 diary, 2 AA batteries, 1 razor 1 Abrading Pad, 2 alcohol wipes, 7 electrodes
2-100-0091	Prep Kit - Midmark IQholter®, 7-lead (case) Each kit contains: 1 diary, 2 AA batteries, 1 razor 1 Abrading Pad, 2 alcohol wipes, 7 electrodes (case - 30 kits)
90-30-0205	Lead Wire, 5-lead, Midmark IQholter® (Blue)
90-30-0206	Lead Wire, 7-lead, Midmark IQholter® (Black)
2-100-0060	Midmark IQholter® Electrodes (box of 50)
2-100-0050	Disposable IQholter® Recorder Pouches (pack of 25)
2-100-0051	Reusable IQholter® Recorder Pouch with Belt



Midmark IQcart® Diagnostic Equipment Cart

(Some assembly required)

- 1 - Storage Basket
- 1 - IQecg® Arm
- 1 - IQspiro® Holster
- 1 - Set of Locking Disks

Model

3-004-1000
3-004-2008
3-004-6028

Description

IQcart® with Arm
IQvitals® Equipment Pole Mount
IQcart® Basket (8" x 6" x 8")



*shown with optional large
Utility Basket

Midmark Care Exchange® 6231 Procedure Workstation with Equipment Pole and Rear Basket

- 29" - 47" Height adjustment
- 20" W x 20" D Worksurface
- 20" W x 18" D Base Dimension
- Weight Capacity: 30 lb
- Unit weight: 37 lb
- Shipping weight: 62 lb
- Equipment pole included
- Rear basket included
- Cuff holder included
- Left and Right Spiro holsters included
- Left and Right Cord Disc's included

Model

6231-001-xxx¹

Description

Procedure Workstation

¹ - See color code chart on page 71

Accessories

Plant Installed	Field Installed	
62XX-501-38	9A489001	Scale Holder
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-04	9A439001	Utility Basket
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-21	9A466001	Large Utility Basket
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W
n/a	3-004-2008	IQvitals® Equipment Pole Mount

**For technical and support information,
call Midmark Technical Service at: 1-800-624-8950, option 2.**



6201-002-xxx shown with optional
coiled power cord w/ (3) outlets

Midmark Care Exchange® 6201 Traditional Style Laptop Workstation

- 29" to 47" height adjustment
- Keyed laptop security system
- 20" D x 20" W worksurface
- 9" D x 20" W secondary shelf
- 9" D x 8" W mouse surface
- 20" D x 18" W base dimension
- Weight capacity: 30 lb
- Unit weight: 44 lb
- Shipping weight: 69 lb

Model

6201-001-xxx¹
6201-002-xxx¹

¹ - See color code chart on page 71

Description

Traditional Style Laptop Workstation
Traditional Laptop Workstation w/Cord Tray

Accessories

Plant Installed	Field Installed	
**62XX-501-02	9A437001	Cable Management Holder
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-07	9A468001	Left Hand Swing-a-Mouse
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
**62XX-501-13	9A451001	Lockable Drawer
n/a	9A460001	Single Flat Panel Monitor Arm
62XX-501-21	9A466001	Large Utility Basket
**62XX-501-22	9A467001	Rear Basket
62XX-501-23	n/a	Med Drawer w/ Lock
62XX-501-24	n/a	Med Drawer No Lock
62XX-501-25	n/a	Med Drawer Chart Pocket
62XX-501-29	9A480001	Cordless Bar Code Holder
*62XX-501-31	9A482001	Cord Tray Rear Basket
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-34	9A485001	Upper Work Surface Extension Brackets
62XX-501-39	9A490001	Cart Handle
*62XX-501-40	9A493001	Coiled Power Cord w/(3) Outlets
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	IQvitals® Equipment Pole Mount
**62XX-501-45	9A594001	Coiled Power Cord
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W

*Only available on 6201-002-XXX models

**Only available on 6201-001-XXX models



Midmark Care Exchange® 6202 Flat Panel PC Workstation

- 29" to 47" height adjustment
- Slim-line PC holder*
- Monitor arm capacity:
10 lb max -001; 20 lb max -002
- 20" D x 20" W worksurface
- 9" D x 8" W mouse surface
- 20" D x 18" W base dimension
- Weight capacity: 30 lb
- Unit weight: 51 lb
- Shipping weight: 76 lb
- VESA mount 75 mm and 100 mm

Model

6202-001-xxx¹
6202-002-xxx¹

¹ - See color code chart on page 71

Description

Flat Panel PC Workstation
HD Flat Panel Workstation

Accessories

Plant Installed	Field Installed	
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-07	9A468001	Left Hand Swing-a-Mouse
62XX-501-08*	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-13	9A451001	Lockable Drawer
62XX-501-19	9A452001	Cable Management Holder
62XX-501-22	9A467001	Rear Basket
62XX-501-26	n/a	Med Drawer w/ Lock
62XX-501-27	n/a	Med Drawer No Lock
62XX-501-28	n/a	Med Drawer Chart Pocket
62XX-501-29	9A480001	Cordless Bar Code Holder
62XX-501-30	9A481001	Monitor Arm TC CPU Holder
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-39	9A490001	Cart Handle
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	IQvitals® Equipment Pole Mount
62XX-501-45	9A594001	Coiled Power Cord
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W

*Standard on 6202-001-xxx only



6203-002-xxx shown with optional coiled power cord w/ (3) outlets

Midmark Care Exchange® 6203 Flat Panel EMR Transition Workstation

- 29" to 47" height adjustment
- Keyed laptop security system
- 23" D x 20" W worksurface
- 12" D x 20" W secondary shelf
- 9" D x 8" W mouse surface
- Monitor arm capacity: 10 lb max
- 20" D x 18" W base dimension
- Weight capacity: 30 lb
- Unit weight: 51 lb
- Shipping weight: 76 lb
- VESA mount 75mm and 100mm

Model

6203-001-xxx¹
6203-002-xxx¹

Description

Flat Panel EMR Transition Workstation
Flat Panel EMR Workstation w/Cord Tray

¹ - See color code chart on page 71

Accessories

Plant Installed	Field Installed	
**62XX-501-02	9A437001	Cable Management Holder
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-07	9A468001	Left Hand Swing-a-Mouse
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
**62XX-501-13	9A451001	Lockable Drawer
62XX-501-21	9A466001	Large Utility Basket
62XX-501-23	n/a	Med Drawer w/ Lock
62XX-501-24	n/a	Med Drawer No Lock
62XX-501-25	n/a	Med Drawer Chart Pocket
62XX-501-29	9A480001	Cordless Bar Code Holder
62XX-501-30	9A481001	Monitor Arm TC CPU Holder
*62XX-501-31	9A482001	Cord Tray Rear Basket
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-34	9A485001	Upper Work Surface Extension Brackets
62XX-501-39	9A490001	Cart Handle
*62XX-501-40	9A493001	Coiled Power Cord w/(3) Outlets
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	IQvitals® Equipment Pole Mount
**62XX-501-45	9A594001	Coiled Power Cord
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W

*Only available on 6203-002-XXX models

**Only available on 6203-001-XXX models



Midmark Care Exchange® 6204 Compact Rectangular Workstation

- 29" to 47" height adjustment
- 15" D x 26" W worksurface
- 20" D x 18" W base dimension
- Weight capacity: 10 lb
- Unit weight: 31 lb
- Shipping weight: 56 lb

Model

6204-001-xxx¹

¹ - See color code chart on page 71

Description

Compact Rectangular Workstation

Accessories

Plant Installed	Field Installed	
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-19	9A452001	Cable Management Holder
62XX-501-21	9A466001	Large Utility Basket
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	IQvitals® Equipment Pole Mount
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W



6205-001-xxx



6205-002-xxx

Midmark Care Exchange® 6205 Compact Kidney Style Workstation

- 29" to 47" height adjustment
- 18" D x 28" W worksurface
- 20" D x 18" W base dimension
- Monitor arm capacity: 10lb max on -002
- Weight capacity: 10 lb
- Unit weight: 34 lb
- Shipping weight: 59 lb

Model

6205-001-xxx¹
6205-002-xxx¹

Description

Compact Kidney Workstation
Flat Panel Compact Kidney Workstation

¹ - See color code chart on page 71

Accessories

Plant Installed	Field Installed	
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-19	9A452001	Cable Management Holder
62XX-501-21	9A466001	Large Utility Basket
*62XX-501-30	9A481001	Monitor Arm CPU Holder
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	IQvitals® Equipment Pole Mount
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W

*Only available on 6205-002-XXX models



Midmark Care Exchange® 6206 Kidney Style Laptop Workstation with Security System

- 29" to 47" height adjustment
- Keyed laptop security cable system
- 22" D x 32" W worksurface
- 20" D x 18" W base dimension
- Weight capacity: 30 lb
- Unit weight: 44 lb
- Shipping weight: 69 lb

Model

6206-001-xxx¹

¹ - See color code chart on page 71

Description

Kidney Style Laptop Workstation with Security System

Accessories

Plant Installed	Field Installed	
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-20	9A452001	Cable Management Holder
62XX-501-21	9A466001	Large Utility Basket
62XX-501-29	9A480001	Cordless Bar Code Holder
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	IQvitals® Equipment Pole Mount
62XX-501-45	9A594001	Coiled Power Cord
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W



Midmark Care Exchange® 6207

Basic Kidney Style Laptop Workstation

- 29" to 47" height adjustment
- 22" D x 32" W worksurface
- 20" D x 18" W base dimension
- Weight capacity: 10 lb
- Unit weight: 37 lb
- Shipping weight: 62 lb

Model

6207-001-xxx¹

¹ - See color code chart on page 71

Description

Kidney Style Laptop Workstation

Accessories

Plant Installed	Field Installed	
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-21	9A466001	Large Utility Basket
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W



Midmark Care Exchange® 6208 Flat Panel Kidney Style Workstation

- 29" to 47" height adjustment
- 22" D x 32" W worksurface
- 20" D x 18" W base dimension
- Monitor arm capacity:
10 lb max -001; 20 lb max -002
- Weight capacity: 30 lb
- Unit weight: 51 lb
- Shipping weight: 76 lb

Model

6208-001-xxx¹
6208-002-xxx¹

¹ - See color code chart on page 71

Description

Flat Panel Kidney Style Workstation
HD Flat Panel Kidney Style Workstation

Accessories

Plant Installed	Field Installed	
62XX-501-04	9A439001	Utility Basket
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-06	9A441001	Bar Code Holder
62XX-501-08	9A443001	Slim Line CPU Holder
62XX-501-09	9A444001	Wide CPU Holder
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-19	9A452001	Cable Management Holder
62XX-501-21	9A466001	Large Utility Basket
62XX-501-29	9A480001	Cordless Bar Code Holder
62XX-501-30	9A481001	Monitor Arm CPU Holder
62XX-501-32	9A483001	Thin Client CPU Holder
62XX-501-41	9A494001	Cart Counterweight
62XX-503-03	n/a	Heavy Duty Casters
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	IQvitals® Equipment Pole Mount
62XX-501-45	9A594001	Coiled Power Cord
62XX-501-38	9A489001	Scale Holder
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W



*shown with optional large Utility Basket

Midmark Care Exchange® 6231 Procedure Workstation with Equipment Pole and Rear Basket

- 29"- 47" Height adjustment
- 20" W x 20" D Worksurface
- 20" W x 18" D Base Dimension
- Weight Capacity: 30 lb
- Unit weight: 37 lb
- Shipping weight: 62 lb
- Equipment pole included
- Rear basket included
- Cuff holder included
- Left and Right Spiro holsters included
- Left and Right Cord Disc's included

Model

6231-001-xxx¹

Description

Procedure Workstation

¹ - See color code chart on page 71

Accessories

Plant Installed	Field Installed	
62XX-501-38	9A489001	Scale Holder
62XX-501-05	9A440001	Multi-Purpose Shelf 20" W
62XX-501-04	9A439001	Utility Basket
62XX-501-12	9A450001	Sharps Container & Bracket
62XX-501-21	9A466001	Large Utility Basket
62XX-501-49	9A596001	Multi-Purpose Shelf Narrow 18" W
n/a	3-004-2008	IQvitals® Equipment Pole Mount



Midmark Care Exchange® 6261 AC Powered Traditional Style Laptop Workstation

- Up to 12 hours of run time between charges
- 29" to 47" height adjustment
- Keyed laptop security system
- 20" D x 20" W worksurface
- 9" D x 20" W secondary shelf
- Coiled power cord
- 9" D x 8" W mouse surface
- 20" D x 18" W base dimension
- Weight capacity: 30 lb
- Unit weight: 90 lb for LiFe option and 127 lb for SLA option
- Shipping weight: 115 lb for LiFe option and 143 for SLA option

Model

6261-001-XXX¹

6261-002-XXX¹

¹ - See color code chart on page 71

Description

AC Traditional Laptop Cart w/LiFe Battery

AC Traditional Laptop Cart w/SLA Battery

Accessories

Plant Installed	Field Installed	
62XX-501-06	9A441001	Bar Code Holder
62XX-501-07	9A468001	Left Hand Swing-a-Mouse
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-23	n/a	Med Drawer w/ Lock
62XX-501-24	n/a	Med Drawer No Lock
62XX-501-25	n/a	Med Drawer Chart Pocket
62XX-501-29	9A480001	Cordless Bar Code Holder
62XX-501-31	9A482001	Cord Tray Rear Basket
62XX-501-34	9A485001	Upper Work Surface Extension Brackets
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	IQvitals® Equipment Pole Mount
62XX-501-38	9A489001	Scale Holder



Midmark Care Exchange® 6263 AC Powered Flat Panel PC Workstation

- Up to 12 hours of run time between charges
- 29" to 47" height adjustment
- Monitor arm: 10 lb max -001 & -002;
20 lb max -003 & -004
- 20" D x 20" W worksurface
- Coiled power cord
- 9" D x 8" W mouse surface
- 20" D x 18" W base dimension
- Weight capacity: 30 lb
- Unit weight: 93 lb for LiFe option
and 130 lb for SLA option
- Shipping weight: 118 lb for
LiFe option and 146 for SLA option
- VESA mount 75 mm and 100 mm

Model

6263-001-XXX¹
6263-002-XXX¹
6263-003-XXX¹
6263-004-XXX¹

Description

AC Flat Panel PC Cart w/LiFe Battery
AC Flat Panel PC Cart w/SLA Battery
AC HD Flat Panel PC Cart w/LiFe Battery
AC HD Flat Panel PC Cart w/SLA Battery

¹ - See color code chart on page 71

Accessories

Plant Installed	Field Installed	
62XX-501-06	9A441001	Bar Code Holder
62XX-501-07	9A468001	Left Hand Swing-a-Mouse
62XX-501-46	9A595001	Equipment Pole / ECG Holder (Rake Style)
62XX-501-23	n/a	Med Drawer w/ Lock
62XX-501-24	n/a	Med Drawer No Lock
62XX-501-25	n/a	Med Drawer Chart Pocket
62XX-501-29	9A480001	Cordless Bar Code Holder
62XX-501-30	9A481001	Monitor Arm TC CPU Holder
62XX-501-31	9A482001	Cord Tray Rear Basket
62XX-501-34	9A485001	Upper Work Surface Extension Brackets
n/a	3-004-6030	Spiro Holster - for Equipment Pole
n/a	3-004-2008	IQvitals® Equipment Pole Mount
62XX-501-38	9A489001	Scale Holder



Midmark Care Exchange® 6281 Laptop Wall Mounted Workstation

- 22" Height adjustment
- 19.5" W x 14" D Worksurface
- Slide-out Mouse Surface
- Up to 45" Reach off wall
- Weight Capacity: 0-8 lb
- Folds to 9" flat to wall
- Unit weight: 33 lb
- Shipping weight: 40 lb

Model	Description
6281-001	Laptop Wall Mounted Workstation

Accessories

Plant Installed	Field Installed	
n/a	9A588001	20" Extension Arm
n/a	9A591001	Track Mount Barcode Holder
n/a	9A590001	Tamper Proof CPU Holder
n/a	9A589001	Adjustable CPU Holder



Midmark Care Exchange® 6282 Standard Duty Wall Mounted Workstation

- 22" Worksurface Height adjustment
- 25" W x 10" D Worksurface
- Max monitor capacity is 20 lb
- 9.5" Independent Monitor Height Adjustment
- Independent 180 (Degree) Monitor rotation and 0-45 (Degree) tilt
- VESA mount 75 mm and 100 mm
- Folds to 9.47" flat to wall
- Weight Capacity: 5-20 lb
- Product weight: 43 lb
- Shipping weight: 50 lb

Model	Description
6282-001	Standard Duty Wall Mounted Workstation

Accessories

Plant Installed	Field Installed	
n/a	9A593001	Dual Monitor Arm
n/a	9A588001	20" Extension Arm
n/a	9A592001	VESA Mount Barcode Holder
n/a	9A591001	Track Mount Barcode Holder
n/a	9A590001	Tamper Proof CPU Holder
n/a	9A589001	Adjustable CPU Holder



Midmark Care Exchange® 6283

Heavy Duty Wall Mounted Workstation

- 22" Worksurface Height adjustment
- 25" W x 10" D Worksurface
- Max monitor capacity is 30 lb
- 9.5" Independent Monitor Height Adjustment
- Independent 180 (Degree) Monitor rotation and 0-45 (Degree) tilt
- VESA mount 75 mm and 100 mm
- Folds to 9.47" flat to wall
- Weight Capacity: 20-35 lb
- Product weight: 43 lb
- Shipping weight: 50 lb

Model

6283-001

Description

Heavy Duty Wall Mounted Workstation

Accessories

Plant Installed	Field Installed	
n/a	9A593001	Dual Monitor Arm
n/a	9A588001	20" Extension Arm
n/a	9A592001	VESA Mount Barcode Holder
n/a	9A591001	Track Mount Barcode Holder
n/a	9A590001	Tamper Proof CPU Holder
n/a	9A589001	Adjustable CPU Holder



Midmark Care Exchange® 6271

Wall Mount Adj. Flip-Up Keyboard/Adj. Flat Screen Monitor

- 29" to 47" height adjustment
- 9" D x 25" W worksurface
- Max monitor weight: 20 lb
- 14" monitor height adjustment
- Worksurface depth adjustable from 14"-18"
- Folds to 5.5" flat to wall
- VESA mount 75 mm and 100 mm
- Product weight: 30 lb
- Shipping weight: 45 lb

Model

6271-001-XXX¹

¹ - See color code chart on page 71

² - Shown with optional shelf

Description

Wall Mount Adj. Flip-Up Keyboard/
Adj. Flat Screen Monitor

Accessories

Plant Installed	Field Installed	
n/a	9A439001	Utility Basket
62XX-501-06	9A441001	Bar Code Holder
n/a	9A443001	Slim Line CPU Holder
n/a	9A444001	Wide CPU Holder
n/a	9A483001	Thin Client CPU Holder
62XX-501-33	9A484001 ¹	Multi-Purpose Shelf
62XX-501-36	9A487001	CPU Holder



Midmark Care Exchange® 6272

Wall Mount Adj. Flat Screen w/flip-up keyboard

- Fixed height
- 9" D x 25" W worksurface
- Max monitor weight is 20 lb
- 14" Monitor height adjustment
- Worksurface depth adjustable from 14"-18"
- Folds to 5.5" flat to wall
- VESA mount 75 mm and 100 mm
- Product weight: 15 lb
- Shipping weight: 23 lb

Model

6272-001-XXX¹

Description

Wall Mount Adj. Screen, Flip-Up Keyboard
Monitor Mount

¹ - See color code chart on page 71

Accessories

Plant Installed	Field Installed	
n/a	9A439001	Utility Basket
n/a	9A443001	Slim Line CPU Holder
n/a	9A483001	Thin Client CPU Holder
62XX-501-33	9A484001 ¹	Multi-Purpose Shelf
62XX-501-36	9A487001	CPU Holder



Midmark Care Exchange® 6273

Wall Mount Monitor Bracket

- VESA mount 75 mm x 100 mm
- Product weight: 1 lb
- Shipping weight: 2 lb
- Monitor weight capacity is 20 lb

Model

6273-001

Description

Wall Mount Bracket



Midmark Care Exchange® 6274

Wall Mount Flip-up Keyboard

- Worksurface depth adjustable from 14"-18" depth
- Folds to 5.5" from wall
- Product weight: 8 lb
- Shipping weight: 11 lb

Model

6274-001-XXX¹

¹ - See color code chart on page 71

Description

Wall Mount Flip-Up Keyboard

Accessories

Plant Installed	Field Installed	
n/a	9A439001	Utility Basket
n/a	9A443001	Slim Line CPU Holder
n/a	9A483001	Thin Client CPU Holder
n/a	9A487001	Wall Column CPU Holder



Midmark Care Exchange® 6275

Wall Mount Adj. Flat Screen Monitor Mount

- 14" Adj. height
- Rotate and tilt
- VESA mount 75 mm and 100 mm
- Monitor weight capacity: 20 lb
- Product weight: 7 lb
- Shipping weight: 10 lb

Model

6275-001

Description

Wall Mount Adj. Flat Screen Monitor
Monitor Mount

Accessories

Plant Installed	Field Installed	
n/a	9A439001	Utility Basket
n/a	9A443001	Slim Line CPU Holder
n/a	9A483001	Thin Client CPU Holder
n/a	9A487001	Wall Column CPU Holder



Midmark Care Exchange® 6276-001

Wall Mount Adj. height w/flip-up keyboard

- 29"- 47" Height adjustment
- 13"x 25" Work surface
- Weight capacity: 20 lb
- Folds to 5.5" to wall
- Product weight: 25 lb
- Shipping weight: 33 lb

Model

6276-001-XXX¹

¹ - See color code chart on page 71

Description

Wall Mount Adj. Flip-Up Keyboard

6276 Accessories

Plant Installed	Field Installed	
n/a	9A439001	Utility Basket
62XX-501-06	9A441001	Bar Code Holder
n/a	9A443001	Slim Line CPU Holder
n/a	9A444001	Wide CPU Holder
n/a	9A483001	Thin Client CPU Holder

Synthesis® Casework Collection by Midmark: Limited Warranty

SCOPE OF WARRANTY

Midmark Corporation warrants to the original purchaser its Midmark Modular Casework products (Synthesis® Casework Collection) manufactured by Midmark to be free from defects in material and workmanship under normal interior use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of the parts or the products the defects of which are reported to Midmark within the applicable warranty period and which upon examination by Midmark prove to be defective.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows:

- A. Five years for all Midmark brand casework products and components (except for the shorter periods in B and C), including door and drawer fronts and casters
- B. Three years for task lights
- C. One year for countertop surfaces and accessories

EXCLUSIONS

This warranty does not cover sinks, faucets and plumbing accessories. Any warranties on these items are extended directly by the manufacturer of these items to the original purchaser. Information on these manufacturers' warranties will be enclosed with the applicable products. Also, Midmark will furnish copies of any warranties extended by any such manufacturers upon request.

This warranty does not cover and Midmark shall not be liable for: (1) repairs and replacements required because of misuse, abuse, negligence, alteration, accident, freight damage or tampering; (2) matching of color, grain or texture except to commercially acceptable standards; (3) changes in color caused by natural or artificial light; (4) products which are not installed, used or properly cleaned as required in Midmark's written Installation Manual; (5) specially manufactured products; (6) products considered to be of a consumable nature such as light bulbs and surge suppression products; (7) accessories or parts not manufactured by Midmark; (8) charges by anyone (including Midmark's authorized dealers) for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which is not expressly authorized in writing in advance by Midmark.

EXCLUSIVE REMEDY

Midmarks' only obligation under this warranty is the repair or replacement of defective parts. Midmark shall not be liable for any direct, special, indirect, incidental, exemplary or consequential damages or delay.

NO AUTHORIZATION

No person or firm is authorized to create for Midmark any other obligation or liability in connection with the products.

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

*Warranty information included in product manuals

Midmark Synthesis® Casework Color-Specific Part Number Explanation and Standard Finish Options

E1RX . 803 . 362

Specific Model Number & Options

Casework Finish

Countertop Color

Example: E1 36" H Exam Room Package, Right Hand Sink, No Options, Henna, Concrete Stone

N2A . 803

Specific Model Number & Options

Casework Finish

Example: Dressing Nook, LH Door, Keyed Locks, Henna



Flare Handle
Serenity Panel

Frame / Panel Color:

Solids:

216 Pebble
490 Fawn
649 Whisper
791 Frost

Metallics:

804 Venus Silver (M)
805 Pearl Essence (M)
806 Earthen Bronze (M)
807 Flat Iron (M)

Others:

793 Pebble PVC-Free

Woodgrains:

459 Acorn
487 Earth
794 Dune
795 Flax
796 Nest
797 Hewn
798 Path
799 Sediment
800 Radiance
801 Timber
802 Storm
803 Henna

Laminate Countertop Color:

362 Concrete Stone
443 Kalahari Topaz
364 Graphite Nebula
360 Evening Tigris
453 Milano Quartz
290 Green Tigris
354 Bronze Legacy
274 Pebble Grey
711 Colorado Slate



Exam Room Casework Package E1

Exam Base: 36" H (39" to top of backsplash) x 48" W x 18" D
(2) Doors w/Waste Drop, (2) 4" Drawers and (2) Doors

Exam Overhead: 48" W x 24" H x 13" D w/3 position Glove,
Towel & Cup Dispenser

- 49" W x 19" D Postform laminate top
- Stainless steel drop-in sink
- Gooseneck facet w/ 4" wrist-blade handles
- Waste can
- Serenity panel style w/flare handles
- Shipping Weight: 327 lb

Product Information

Step 1: Configuration

E1

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Casework Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E1LX.803.362

*See page 71 for Casework Finish & Countertop Color



Exam Room Casework Package E2

Exam Base: 36" H (39" to top of backsplash) x 48" W x 18" D
(2) Doors w/Waste Drop, (3) 4" Drawers and (1) 8" Drawer

Exam Overhead: 48" W x 24" H x 13" D w/3 position Glove,
Towel & Cup Dispenser

- 49" W x 19" D Postform laminate top
- Stainless steel drop-in sink
- Gooseneck facet w/ 4" wrist-blade handles
- Waste can
- Serenity panel style w/flare handles
- Shipping Weight: 331 lb

Product Information

Step 1: Configuration

E2

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Casework Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E2LX.803.362

*See page 71 for Casework Finish & Countertop Color



Exam Room Casework Package E3

Overall Exam Base Dimensions: 36" H x 54" W x 18" D w/Support Leg

Exam Base Cabinet: 36" H (39" to top of backsplash) x 24" W x 18" D
(2) Doors w/Waste Drop

Exam Cart: 33" H x 24" W x 18" D (3) 4" Drawers, (1) 8" Drawer
w/Kydex Contour top and 60mm casters

Exam Overhead: 24" H x 48" W x 13" D w/(3) position Glove,
Towel & Cup Dispenser



- 55" W x 19" D Postform laminate top
- Stainless steel drop-in sink
- Gooseneck facet w/ 4" wrist-blade handles
- Waste can
- Serenity panel style w/flare handles
- Shipping Weight: 437 lb

Product Information

Step 1: Configuration

E3

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Casework Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E3LX.803.362

*See page 71 for Casework Finish & Countertop Color



Exam Room Casework Package E4

Overall Exam Base Dimensions: 36" H (39" to top of backsplash)
x 48" W x 18" D

Exam Base Cabinet: 36" H x 24" W x 18" D (2) 4" Drawers & (2) Doors
Wall Mounted Sink: 24" W X 18" D

Exam Overhead: 24" H & 48" W x 13" D w/(3) position Glove,
Towel & Cup Dispenser

- 49" W x 19" D Postform laminate top
- Stainless steel drop-in wall mounted sink
- Gooseneck facet w/ 4" wrist-blade handles
- Serenity panel style w/flare handles
- Shipping Weight: 307 lb

Product Information

Step 1: Configuration

E4

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Casework Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E4LX.791.362

*See page 71 for Casework Finish & Countertop Color



Exam Room Casework Package E5

Overall Exam Base Dimensions: 36" H (39" to top of backsplash)
x 54" W x 18" D w/Support Leg

Exam Cart: 33" H x 24" W x 18" D Cart (3) 4" Drawers, (1) 8" Drawer
w/Kydex Contour Top and 60mm Casters
Wall Mounted Sink: 24" W X 18" D

Exam Overhead: 24" H & 48" W x 13" D w/(3) position Glove,
Towel & Cup Dispenser

- 55" W x 19" D Postform laminate top
- Stainless steel drop-in wall mounted sink
- Gooseneck facet w/ 4" wrist-blade handles
- Serenity panel style w/flare handles
- Shipping Weight: 401 lb

Product Information

Step 1: Configuration

E5

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Casework Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E5LX.803.362

*See page 71 for Casework Finish & Countertop Color



Exam Room ADA Casework Package E6

Overall Exam Base Dimensions: 34" H ADA (37" to top of backsplash)
x 48" W x 18" D

Exam Base Cabinet: 34" H ADA x 18" W x 18" D (2) 4" Drawers
& (1) 8" Drawer

ADA Wall Mounted Sink: 30" W X 18" D

Exam Overhead: 24" H x 48" W x 13" D w/(3) position Glove,
Towel & Cup Dispenser

- 49" W x 19" D Postform laminate top
- Stainless steel drop-in wall mounted sink
- Gooseneck facet w/ 4" wrist-blade handles
- Serenity panel style w/flare handles
- Shipping Weight: 294 lb

*Customer is responsible for ensuring product meets Federal, state and local code.

Product Information

Step 1: Configuration

E6

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Casework Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E6LX.803.362

*See page 71 for Casework Finish & Countertop Color



Exam Room ADA Casework Package E7

Overall Exam Base Dimensions: 34" H ADA (37" to top of backsplash)
x 60" W x 18" D w/Support Leg

Exam Cart: 31" H x 24" W x 18" D Cart (2) 4" Drawers, (1) 8" Drawer
w/Kydex® Contour Top and 60 mm Casters

Wall Mounted Sink: 30" W X 18" D

Exam Overhead: 24" H x 48" W x 13" D w/(3) position Glove,
Towel & Cup Dispenser

- 61" W x 19" D Postform laminate top
- Stainless steel drop-in wall mounted sink
- Gooseneck facet w/ 4" wrist-blade handles
- Serenity panel style w/flare handles
- Shipping Weight: 294 lb

*Customer is responsible for ensuring product meets Federal, state and local code.

Product Information

Step 1: Configuration

E7

Step 2: Sink Location

L Left Hand

R Right Hand

Step 3: Option Package

X None

A Keyed Lock

Step 4: Casework Finish*

XXX Solid Color Cabinet Finish
(Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free
Cabinet Finish

Step 5: Countertop Color*

XXX Countertop Color

Configuration Example: E7LX.803.362

*See page 71 for Casework Finish & Countertop Color



* Shown with Kydex®
contour top

Mobile Treatment Cabinet M2

33" H x 27" W x 18" D (2) 4" Drawers & (2) Doors

- Serenity panel style w/flare handles
- Shipping Weight: 155 lb

Product Information

Step 1: Configuration

M2

Step 2. Depth

1 18" Deep

Step 3. Option Package

- A** 60mm Casters
- B** 4" Quiet Casters
- C** 60mm Casters & Keyed Locks
- D** 4" Quiet Casters & Keyed Locks

Step 4. Casework Finish*

- XXX** Solid Color Cabinet Finish (Excludes PVC-Free)
- XXX** Woodgrain, Metallic & PVC-Free Cabinet Finish

Step 5. Countertop Type

- KCN** Kydex® Contour Top
- SEB** Soft Edge Bumper Top
- SEH** Soft Edge Handle Top (29"W)

Configuration Example: M21A.791.KCN

*See page 71 for Casework Finish



* Shown with soft edge bumper top

Mobile Treatment Cabinet M4

33" H x 27" W x 18" D (3) 4" Drawers & (1) 8" Drawer

- Serenity panel style w/flare handles
- Shipping Weight: 152 lb

Product Information

Step 1: Configuration

M4

Step 2. Depth

1 18" Deep

Step 3. Option Package

- A** 60mm Casters
- B** 4" Quiet Casters
- C** 60mm Casters & Keyed Locks
- D** 4" Quiet Casters & Keyed Locks

Step 4. Casework Finish*

- XXX** Solid Color Cabinet Finish (Excludes PVC-Free)
- XXX** Woodgrain, Metallic & PVC-Free Cabinet Finish

Step 5. Countertop Type

- KCN** Kydex® Contour Top
- SEB** Soft Edge Bumper Top
- SEH** Soft Edge Handle Top (29"W)

Configuration Example: M41A.791.SEB

*See page 71 for Casework Finish



Mobile Treatment Cabinet M5

33" H x 27" W x 18" D (5) 4" Drawers

- Serenity panel style w/flare handles
- Shipping Weight: 159 lb

* Shown with soft edge
handle top

Product Information

Step 1: Configuration

M5

Step 2. Depth

1 18" Deep

Step 3. Option Package

- A** 60mm Casters
- B** 4" Quiet Casters
- C** 60mm Casters & Keyed Locks
- D** 4" Quiet Casters & Keyed Locks

Step 4. Casework Finish*

- XXX** Solid Color Cabinet Finish (Excludes PVC-Free)
- XXX** Woodgrain, Metallic & PVC-Free Cabinet Finish

Step 5. Countertop Type

- KCN** Kydex® Contour Top
- SEB** Soft Edge Bumper Top
- SEH** Soft Edge Handle Top (29"W)

Configuration Example: M51A.791.SEH

*See page 71 for Casework Finish

Glove, Towel & Cup Dispenser



* D2 shown

- 18" W x 8" D x 24" H
- Shipping weight: 37 lb
- Serenity panel style w/flare handle



* D4 shown

Product Information

Step 1. Configuration

D1 Glove/Towel/Cup Dispenser - (2) Opening RH Door

D2 Glove/Towel/Cup Dispenser - (2) Opening LH Door

D3 Glove/Towel/Cup Dispenser - (3) Opening RH Door

D4 Glove/Towel/Cup Dispenser - (3) Opening LH Door

Step 2. Option Package

X No Lock

Step 3. Casework Finish*

XXX Solid Color Cabinet Finish (Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free Cabinet Finish

Configuration Example: D2X.803

*See page 71 for Casework Finish



Dressing Nook

- Left-hand door pictured
- 24" W x 7" D x 63" H
- Interior mirror, garment hooks and shoe storage included
- Serenity panel style w/flare handle
- Shipping weight: 138 lb



Product Information

Step 1. Configuration

N1 Dressing Nook - RH Door

N2 Dressing Nook - LH Door

Step 2. Option Package

X No Lock

A Keyed Locks

Step 3. Casework Finish*

XXX Solid Color Cabinet Finish (Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free Cabinet Finish

Configuration Example: N2X.803

*See page 71 for Casework Finish



Waste Receptacle

- 12" W x 12" D x 24" H
- Capacity: 33 liquid qts
- Serenity panel style
- Shipping weight: 35 lb

Product Information

Step 1. Configuration

W1 Waste Receptacle

Step 2. Option Package

X No Lock

Step 3. Casework Finish*

XXX Solid Color Cabinet Finish (Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free Cabinet Finish

Configuration Example: W1X.803

*See page 71 for Casework Finish



* I2 shown

Instrument Panel

- Instruments not included
- Shipping weight: 17 lb

Product Information

Step 1. Configuration

I1 Instrument Panel - 30" x 20"

I2 Instrument Panel - 34" x 20"

Step 2. Option Package

X No Lock

Step 3. Casework Finish*

XXX Solid Color Cabinet Finish (Excludes PVC-Free)

XXX Woodgrain, Metallic & PVC-Free Cabinet Finish

Configuration Example: I2X.803

*See page 71 for Casework Finish



Ritter Waste Receptacle

- With bag retaining ring
- 14" W x 16" D x 25" H
- Shipping weight: 17 lb
- Capacity: 43 liquid qts.
- Color: Pebble

Model

261-001

Description

Waste receptacle



Ritter M11 UltraClave®

Automatic Sterilizer

- Automatic Door
- Chamber size: 11" x 18"
- Meets UL 61010-1, 2nd Edition, IEC 61010-2-040, 1st Edition
- Meets CAN/CSA C22.2 No. 61010-1, 2nd Edition
- Meets CSA C22.2 No. 61010-2-040-07 Part 2-040, 1st Edition
- ASME Certified
- Sterilizer weight: 99 lb (44.9 kg)
- Shipping weight: 131 lb (59.4 kg)

Model

M11-022

Description

Ritter M11 UltraClave® - large table-top automatic sterilizer with automatic door



Ritter M9 UltraClave®

Automatic Sterilizer

- Automatic Door
- Chamber size: 9" x 15"
- Meets UL 61010-1, 2nd Edition, IEC 61010-2-040, 1st Edition
- Meets CAN/CSA C22.2 No. 61010-1, 2nd Edition
- Meets CSA C22.2 No. 61010-2-040-07 Part 2-040, 1st Edition
- ASME Certified
- Sterilizer weight: 73 lb (33.1 kg)
- Shipping weight: 81 lb (36.7 kg)

Model

M9-022

Description

Ritter M9 UltraClave® - fully automatic sterilizer with automatic door



Ritter M9D AutoClave® Manual Door

- Chamber size: 9" x 15"
- Meets UL 61010-1, 2nd Edition, IEC 61010-2-040, 1st Edition
- Meets CAN/CSA C22.2 No. 61010-1, 2nd Edition
- Meets CSA C22.2 No. 61010-2-040-07 Part 2-040 1st Edition
- ASME Certified
- Sterilizer weight: 73 lb (33.1 kg)
- Shipping weight: 81 lb (36.7 kg)

Model

M9D-022

Description

Ritter M9D AutoClave® - automatic sterilizer with manual door



Ritter M3 UltraFast® Automatic Sterilizer

- Automatic Door
- Chamber size: 12.1" L x 7.6" W x 1.6" D (0.49 gal.)
- Meets UL 61010-1 2nd edition, IEC 61010-2-040
- Meets CAN/CSA C22.2 No. 60110-1 2nd Edition
- Meets CSA C22.2 No. 61010-2-040
- FCC Part 15, Sub-part B
- ASME Certified
- Sterilizer weight: 71 lb
- Shipping weight: 80 lb

Model

M3-003

Description

Ritter M3 UltraFast® - compact, fast, automatic sterilizer

Sterilizer Accessories

Model	Description
002-0396-00	Speed-Clean (16 oz. single bottle) (C)
002-0396-05	Speed-Clean (12 - 16 oz. bottles)
9A226001	M9/M11 Pouch Rack (D)
9A259001	M9/M11 Impact Dot Matrix Printer (B)
9A307001	M9/M11 Cool Hand Tool (Tray and Cassette Removal Tool) (A)
9A401001	M3 Printer (E)
9A402001	M3 Door Tray (F)
9A404001	M3 Top Cover Protector (G)

(A) Cool Handtool

8.5" L



(B) Impact Dot Matrix Printer

7.5" L x 4.4" W x 2.25" H




(C) SpeedClean

8" H



(D) Pouch Rack



(E) M3 Printer



(F) M3 Door Tray



(G) M3 Top Cover Protector





Single VistaCool™ Direct-to-Drain System

- Designed to be compatible with the Ritter M3 sterilizer
- Tank dimensions w/fittings: 7.6" W x 7.6" L x 13.75" H
(19.3 cm x 19.3 cm x 34.9 cm)
- Tank volume: 1.50 gal (5.70 liters)
- Tank weight w/water and fittings: 15.0 lb (6.8 kg)
- Approximate shipping weight: 7.0 lb (3.2 kg)
- Max coolant water pressure: 100 psi (6.9 bar)

Model

9A586001

Description

VistaCool™ single unit for connecting to one M3 sterilizer



Double VistaCool™ Direct-to-Drain System

- Designed to be compatible with any two current Ritter sterilizers
- Tank dimensions w/fittings: 4.0" W x 13.0" L x 17.0" H
(10.2 cm x 33.1 cm x 43.2 cm)
- Tank volume: 3.00 gal (11.40 liters)
- Tank weight w/water and fittings: 28.0 lb (12.7 kg)
- Approximate shipping weight: 10.0 lb (4.5 kg)
- Max coolant water pressure: 100 psi (6.9 bar)

Model

9A586002*

Description

VistaCool™ double unit for connecting to any two current Ritter sterilizers

*Required adapter kits included. Ritter M9 and M11 only compatible with the double unit.

Accessories

Model

Description

9A586004

Air Gap Kit (needed only if required by local code)

9A586005

Single VistaCool™ adjustable tank stand

9A586006

Double VistaCool™ adjustable tank stand

VistaCool™ is a trademark of Vista Research Group, Ashland, OH



QuickClean™ Ultrasonic Cleaner QC6/QC6R

- Tank Capacity: 6.6 gal/25L
- Tank Size: 19.9" L x 11.8" D x 7.8" H (50.5 x 30 x 20 cm)
- Overall Size: 21.6" L x 16.5" D x 14.3" H (55 x 42 x 36.5 cm)
- Table Top Unit Weight: 29.6 lb (13.44 kg)
- Table Top Unit Shipping Weight: 32.6 lb (14.78 kg)
- Recessed Unit Weight: 32.9 lb (14.92 kg)
- Recessed Unit Shipping Weight: 35.8 lb (16.26 kg)
- Safety basket INCLUDED!
- CSA Listed
- TUV and FCC Listed

Model	Description
QC6-01	QuickClean™ Ultrasonic Cleaner – 6.6 gal Tabletop Unit
QC6R-01	QuickClean™ Ultrasonic Cleaner – 6.6 gal Recessed Unit



QuickClean™ Ultrasonic Cleaner QC3/QC3R

- Tank Capacity: 3.3 gal/12.5L
- Tank Size: 13" L x 11.8" W x 5.9" H (33 x 30 x 15 cm)
- Overall Size: 14.7" L x 16.5" W x 10.4" H (37.5 x 42 x 26.5 cm)
- Table Top Unit Weight: 18.6 lb (8.42 kg)
- Table Top Unit Shipping Weight: 20.9 lb (9.50 kg)
- Recessed Unit Weight: 22.3 lb (10.10 kg)
- Recessed Unit Shipping Weight: 24.6 lb (11.18 kg)
- Safety basket INCLUDED!
- CSA Listed
- TUV and FCC Listed

Model	Description
QC3-01	QuickClean™ Ultrasonic Cleaner – 3.3 gal Tabletop Unit
QC3R-01	QuickClean™ Ultrasonic Cleaner – 3.3 gal Recessed Unit



QuickClean™ Ultrasonic Cleaner QC1

*This specific unit
will be available
for sale on
7/26/2017.*

- Tank Capacity: 1.2 gal/4.5L
- Tank Size: 11.8" L x 5.9" W x 5.9" H (30 x 15 x 15 cm)
- Overall Size: 13.5" L x 10.2" W x 10.4" H (30 x 26 x 26.5 cm)
- Table Top Unit Weight: 12.7 lb (5.76 kg)
- Table Top Unit Shipping Weight: 14.6 lb (6.64 kg)
- Safety basket INCLUDED!
- CSA Listed
- TUV and FCC Listed

Model	Description
QC1-01	QuickClean™ Ultrasonic Cleaner – 1.2 gal Tabletop Unit

Ultrasonic Parts & Accessories

9A612001	QC1 QuickClean™ Beaker Accessory Kit
9A613001	QC3/QC3R QuickClean™ Beaker Accessory Kit
9A614001	QC6/QC6R QuickClean™ Beaker Accessory Kit
9A298001	Midmark Enzymatic Cleaner (32 oz)
9A296001	Midmark General Purpose Cleaner (32 oz)
9A297001	Midmark Tartar and Stain Remover (32 oz)
002-10007-00	QC1 Safety Basket (1 basket comes standard with a QC1 unit)
002-10008-00	QC3/QC3R Safety Basket (1 basket comes standard with QC3/QC3R unit)
002-10009-00	QC6/QC6R Safety Basket (1 basket comes standard with QC6/QC6R unit)

QUESTIONS? Call us.
1-800-MIDMARK

CARB 93120.2 Phase 2 Compliant

Midmark is an ISO 13485 and ISO 9001 Certified Company.
For more information or a demonstration, contact your Midmark
dealer or call: 1-800-MIDMARK Fax: 1-800-365-8631
or visit our website at midmark.com

© 2008 Midmark Corporation
Midmark Corporation, Dayton, OH.
Products subject to improvement changes without notice
Litho in U.S.A. MKT-00105 Rev. B1 (6/17)



Because we care.